


Solvia Rise DR-T

MR Conditional Dual-Chamber Pacemaker



Ordering Information

Model	Connectors	Volume/Weight	Dimensions	Order number
Solvia Rise DR-T	IS-1 (2x)	11 cm ³ /23.2 g	48 x 44 x 6.5 mm	460227

Product Highlights

Closed Loop Stimulation (CLS)

BIOTRONIK Home Monitoring®

MRI Guard 24/7

EarlyCheck

Solvía Rise DR-T

Technical Data

MR Conditional	
ProMRI	For combination of MR Conditional devices, please see the manual "ProMRI MR Conditional device systems"
Closed Loop Stimulation	
CLS mode	DDD-CLS; VI-CLS
Max. CLS rate	80 ... (10) ... 180 bpm
Expert parameters	
CLS response	Very low; Low; Medium; High; Very high
CLS resting rate control	OFF; +10 ... (10) ... +50 bpm
Vp required	Yes ; No
Pacing Parameters	
NBD-NBG code	00E-DDR
Mode	DDD-CLS; VI-CLS; DDDR-ADIR; DDDR; DDIR; VVIR; AAIR; D00; DDD-ADI; DDD; DDI; VVI; AAI; V00; VDDR; VDIR; DDT; VDD; VDI; WT; OFF
Basic rate	30 ... (5) ... 100 ... (10) ... 200 bpm
Night rate	OFF; 30 ... (5) ... 100 bpm
Rate hysteresis	OFF; -5 ... (-5) ... -25 ... (-20) ... -65 bpm
Scan/Repetitive	OFF; ON
Rate fading	OFF; ON
Mode switching	
Intervention rate	OFF; 100 ... (10) ... 250 bpm
Onset criterion	3 ... (1) ... 8 out of 8
Resolution criterion	3 ... (1) ... 8 out of 8
Mode switching expert parameters	
Change of basic rate	OFF; +5 ... (5) ... +30 bpm
Post ModeSw rate	OFF; +5 ... (5) ... +50 bpm
Post ModeSw duration	1 ... (1) ... 30 min
Rate stabilization during mode switching	OFF; ON
2:1 Lock-in protection	OFF; ON
Vp suppression	OFF; ON (programmable via the modes DDDR-ADIR and DDD-ADI)
Pacing suppression	1 ... (1) ... 8 consecutive Vs
Pacing support	1 ... (1) ... 4 out of 8 cycles
Pulse amplitude (A, RV)	0.5 ... (0.25) ... 4.0 ... (0.5) ... 6.0; 7.5 V
Pulse width (A, RV)	0.1 ... (0.1) ... 0.5 ... (0.25) ... 1.5 ms
Capture control (A, RV)	OFF; ATM; ON
Sensing (A)	OFF; AUTO; 0.1 ... (0.1) ... 1.5 ... (0.5) ... 7.5 mV
Sensing (RV)	OFF; AUTO; 0.5 ... (0.5) ... 7.5 mV
AV delay	15...(5)...300 ms, fixed; 40...(5)...350 ms dynamic
AV dynamics	Low; Medium; High; Fixed
Sense compensation	OFF; -5 ... (-5) ... -120 ms
AV hysteresis mode	OFF; Positive; Negative; IRSplus

Arrhythmia Detection and Redection	
AT/AF detection	HAR limit
HAR limit	100 ... (10) ... 250 bpm
AT/AF counter	Detection counter: 36 out of 48; Termination counter: 20 out of 24
HVR detection	HVR limit
HVR limit	4 ... (4) ... 20 ... (5) ... 60 Events
Tachycardia Therapy	
Atr. NIPS	Burst Pacing; Programmed Stimulation
PMT detection/termination	OFF; ON
Leads	
Automatic lead check (A, RV)	ON; OFF
Lead configuration (A, RV)	Unipolar; Bipolar
MRI	
MRI mode	A00; D00; V00; AUTO
Basic rate	Mean rate + 15 bpm; 70 ... (5) ... 100 ... (10) ... 160 bpm
Timing Parameters	
Upper rate	90 ... (10) ... 200 bpm
Atrial upper rate	OFF; 175; 200; 240 bpm
Atrial refractory period	AUTO (AV delay but at least 225 ms)
RV refractory period	200 ... (25) ... 500 ms
PVARP	AUTO; 175 ... (25) ... 600 ms
Blanking after atr. pacing RV	30 ... (5) ... 70 ms
Physical Parameters	
Service time	14 years, 10 months ¹¹ ¹¹ A/RV: 2.5 V/0.4 ms, 60 bpm, 500 Ω; A: pacing: 50 %; RV: pacing: 5 %; Home Monitoring: OFF, QuickCheck: OFF, RF telemetry: OFF, Vp Suppression: ON
Electrically conductive surface	29.9 cm ²
Additional Parameters	
Magnet response	AUTO (10 cycles asynchronous; then basic rate synchronous), asynchronous, synchronous
IEGM recording	17 recordings, at least 16 seconds each Total recording time approx. 5 minutes
Auto initialization	ON
Recording Episodes	
High atrial rate	OFF; ModeSw; AT
High ventricular rate	OFF; ON
For nsT	OFF; ON
Patient trigger	OFF; ON
Tests	
Different tests for	Impedance, Sensing, Pacing threshold, Retrograde conduction, AV optimization
Program Sets	
Programs	Standard program; ProgramConsult; Individual program (1-3, individually program- mable); First interrogated program; Safe program

BIOTRONIK Home Monitoring®

Transmitted Data	AF diagnostics, Heart Failure Monitor diagnostics, Detection and therapy counters, Statistics, Lead measurement values, Battery status, Program parameters, Threshold (A/RV), Sensing amplitude (A/RV), Pacing statistics, Arrhythmia statistics (A/RV)
------------------	---

Message types

Trend message	Triggered automatically once every 24 hours
---------------	---

Programmer Settings


Home Monitoring	OFF; ON (ON by default)
Ongoing atrial episode	OFF; 6 h; 12 h; 18 h

Home Monitoring-supported follow-up

Remote Scheduling	Enable; Disable
HM follow-up intervals/alignment	Individually programmable first date and repetition intervals varying from 20-366 days; Alignment with a specific day of the week; Only working days or no day alignment
EarlyCheck	Automatic first Home Monitoring-supported follow-up 2 hours after implantation detection
Transmitted data	Periodic IEGM; Rate histogram (V); Device settings and statistics
Please refer to the technical manual of the device for further technical information.	

Solvia Rise SR-T

MR Conditional Single-Chamber Pacemaker



Ordering Information

Model	Connectors	Volume/Weight	Dimensions	Order number
Solvia Rise SR-T	IS-1 (1x)	10 cm ³ /20.8 g	48 x 40 x 6.5 mm	460228

Product Highlights

Closed Loop Stimulation (CLS)

BIOTRONIK Home Monitoring®

MRI Guard 24/7

EarlyCheck

Solvia Rise SR-T

Technical Data

MR Conditional	
ProMRI	For combination of MR Conditional devices, please see the manual "ProMRI MR Conditional device systems"
Closed Loop Stimulation	
CLS mode	VVI-CLS
Max. CLS rate	80 ... (10) ... 180 bpm
Expert parameters	
CLS response	Very low; Low; Medium; High; Very high
CLS resting rate control	OFF; +10 ... (10) ... +50 bpm
Vp required	Yes ; No
Pacing Parameters	
NBD-NBG code	00E-WVIR/00E-AAIR
Lead position	A, RV
Mode (when lead position RV)	VVI-CLS; WVIR; VVI; V00; VVT; OFF
Mode (when lead position A)	AAIR; AAI; OFF
Basic rate	30 ... (5) ... 100 ... (10) ... 200 bpm
Night rate	OFF; 30 ... (5) ... 100 bpm
Rate hysteresis	OFF; -5 ... (-5) ... -25 ... (-20) ... -65 bpm
Scan/Repetitive	OFF; ON
Rate fading	OFF; ON
Pulse amplitude	0.5 ... (0.25) ... 4.0 ... (0.5) ... 6.0; 7.5 V
Pulse width	0.1 ... (0.1) ... 0.5 ... (0.25) ... 1.5 ms
Capture control	OFF; ATM; ON
Sensing (A)	OFF; AUTO; 0.1 ... (0.1) ... 1.5 ... (0.5) ... 7.5 mV
Sensing (RV)	AUTO; 0.5 ... (0.5) ... 7.5 mV
Arrhythmia Detection and Redection	
AT/AF detection (when lead position A)	HAR limit
HAR limit	100 ... (10) ... 250 bpm
HVR detection (when lead position RV)	HVR limit
HVR limit	150 ... (5) ... 200 bpm
HVR counter	4 ... (4) ... 20 ... (5) ... 60 Events
Leads	
Automatic lead check	ON; OFF
Lead configuration	Unipolar; Bipolar

MRI	
MRI mode (when lead position RV)	V00; AUTO
MRI mode (when lead position A)	A00; AUTO
Basic rate	70 ... (5) ... 100 ... (10) ... 160 bpm
Timing Parameters	
Atrial refractory period (when lead position A)	300 ... (25) ... 775 ms
Sensing	AUTO; 0.5 ... (0.5) ... 7.5 mV
RV refractory period (when lead position RV)	200 ... (25) ... 500 ms
Physical Parameters	
Service time	15 years, 1 month ¹⁾ ¹⁾ At 2.5 V/0.4 ms, 60 bpm, 500 Ω; pacing: 50 %; Home Monitoring: OFF, RF telemetry: OFF
Electrically conductive surface	29.9 cm ²
Additional Parameters	
Magnet response	AUTO (10 cycles asynchronous; then basic rate synchronous), asynchronous, synchronous
IEGM recording	17 recordings, at least 16 seconds each Total recording time approx. 5 minutes
Auto initialization	ON
Recording Episodes	
High atrial rate (when lead position A)	OFF; AT
High ventricular rate (when lead position RV)	OFF; ON
For nsT (when lead position RV)	OFF; ON
Patient trigger	OFF; ON
Tests	
Different tests for	Impedance, Sensing, Pacing threshold
Program Sets	
Programs	Standard program; ProgramConsult; Individual program (1-3, individually programmable); First interrogated program; Safe program

BIOTRONIK Home Monitoring®

Transmitted Data	Heart Failure Monitor diagnostics, Detection counters, Statistics, Lead measurement values, Battery status, Program parameters, Threshold, Sensing amplitude, Pacing statistics, Arrhythmia statistics
------------------	--

Message types

Trend message	Triggered automatically once every 24 hours
---------------	---

Programmer Settings

Home Monitoring	OFF; ON (ON by default)
-----------------	-------------------------

Home Monitoring-supported follow-up

Remote Scheduling	Enable; Disable
HM follow-up intervals/alignment	Individually programmable first date and repetition intervals varying from 20-366 days; Alignment with a specific day of the week; Only working days or no day alignment
EarlyCheck	Automatic first Home Monitoring-supported follow-up 2 hours after implantation detection
Transmitted data	Periodic IEGM; Rate histogram (V); Device settings and statistics

Please refer to the technical manual of the device for further technical information.