

# AlCath FullCircle

## The circle closes



### Curve Reference Chart

| Model                    |  | Curve Type        |                    |            |            |                              |         |
|--------------------------|--|-------------------|--------------------|------------|------------|------------------------------|---------|
| <b>BIOTRONIK</b>         | <b>AlCath FullCircle/<br/>AlCath Flux FullCircle</b> | Red               | Blue               | Green      | Black      | Cyan                         |         |
|                          | <b>Tip Length T in mm</b>                            | 65                | 75                 | 85         | 95         | 105                          |         |
|                          | <b>Reach R in mm</b>                                 | 48                | 57                 | 65         | 73         | 80                           |         |
|                          | <b>Curve Radius C in mm</b>                          | 28                | 38                 | 43         | 48         | 56                           |         |
| <b>Biosense Webster</b>  | <b>Celsius/ NaviStar<br/>(ThermoCool)</b>            | C (Green)         | D (Blue)           | E (White)  | F (Orange) | G (Gray) H (Brown) J (Black) |         |
| <b>St. Jude Medical</b>  | <b>Livewire TC</b>                                   | Small             | Medium             |            |            | Large                        |         |
|                          | <b>Therapy</b>                                       | Medium            |                    | Ext. R.    | Far R.     | Large                        | X-Large |
|                          | <b>Therapy Cool Path</b>                             | Medium            |                    | L1         | FL         | Large                        | X-Large |
| <b>Boston Scientific</b> | <b>Blazer II (XP)</b>                                | Small (K1)        | Standard           |            |            | Large (K2)                   |         |
|                          | <b>Chilli II</b>                                     |                   | Standard           |            |            | Large (K2)                   |         |
| <b>Medtronic</b>         | <b>RF Enhancr II</b>                                 |                   | Standard           | Medium     |            | Large                        |         |
|                          | <b>Sprinklr</b>                                      |                   |                    | Standard   |            | X-Large                      |         |
| <b>Bard</b>              | <b>Stinger</b>                                       | B (Red) C (Green) | D (Blue) E (White) | F (Orange) |            | G (Gray)                     |         |

| Model                                   | Curve Type | Electrode Material | Order Number | Model                                 | Curve Type | Electrode Material | Order Number |
|---|------------|--------------------|--------------|---------------------------------------|------------|--------------------|--------------|
| <b>AlCath TC G FullCircle (4 mm)</b>    |            |                    |              | <b>AlCath TC FullCircle (4 mm)</b>    |            |                    |              |
| AlCath Red TC G FullCircle              | Red        | Au                 | 364 537      | AlCath Red TC FullCircle              | Red        | Pt/Ir              | 364 536      |
| AlCath Blue TC G FullCircle             | Blue       | Au                 | 364 525      | AlCath Blue TC FullCircle             | Blue       | Pt/Ir              | 359 038      |
| AlCath Green TC G FullCircle            | Green      | Au                 | 364 533      | AlCath Green TC FullCircle            | Green      | Pt/Ir              | 364 532      |
| AlCath Black TC G FullCircle            | Black      | Au                 | 364 522      | AlCath Black TC FullCircle            | Black      | Pt/Ir              | 359 037      |
| AlCath Cyan TC G FullCircle             | Cyan       | Au                 | 364 529      | AlCath Cyan TC FullCircle             | Cyan       | Pt/Ir              | 364 528      |
| <b>AlCath TC LT G FullCircle (8 mm)</b> |            |                    |              | <b>AlCath TC LT FullCircle (8 mm)</b> |            |                    |              |
| AlCath Red TC LT G FullCircle           | Red        | Au                 | 364 539      | AlCath Red TC LT FullCircle           | Red        | Pt/Ir              | 364 538      |
| AlCath Blue TC LT G FullCircle          | Blue       | Au                 | 364 527      | AlCath Blue TC LT FullCircle          | Blue       | Pt/Ir              | 364 526      |
| AlCath Green TC LT G FullCircle         | Green      | Au                 | 364 535      | AlCath Green TC LT FullCircle         | Green      | Pt/Ir              | 364 534      |
| AlCath Black TC LT G FullCircle         | Black      | Au                 | 364 524      | AlCath Black TC LT FullCircle         | Black      | Pt/Ir              | 364 523      |
| AlCath Cyan TC LT G FullCircle          | Cyan       | Au                 | 364 531      | AlCath Cyan TC LT FullCircle          | Cyan       | Pt/Ir              | 364 530      |
| <b>AlCath TC Flux G FullCircle</b>      |            |                    |              | <b>AlCath TC Flux FullCircle</b>      |            |                    |              |
| AlCath Flux Red TC G FullCircle         | Red        | Au                 | 364 550      | AlCath Flux Red TC FullCircle         | Red        | Pt/Ir              | 364 549      |
| AlCath Flux Blue TC G FullCircle        | Blue       | Au                 | 364 542      | AlCath Flux Blue TC FullCircle        | Blue       | Pt/Ir              | 359 216      |
| AlCath Flux Green TC G FullCircle       | Green      | Au                 | 364 547      | AlCath Flux Green TC FullCircle       | Green      | Pt/Ir              | 364 546      |
| AlCath Flux Black TC G FullCircle       | Black      | Au                 | 364 541      | AlCath Flux Black TC FullCircle       | Black      | Pt/Ir              | 359 215      |
| AlCath Flux Cyan TC G FullCircle        | Cyan       | Au                 | 364 545      | AlCath Flux Cyan TC FullCircle        | Cyan       | Pt/Ir              | 364 543      |

| Patient Cable | Length | Ablation System                   | Order Number |
|---------------|--------|-----------------------------------|--------------|
| PK-185        | 2,5 m  | BIOTRONIK AbControl MD5           | 365 954      |
| PK-111        | 2,0 m  | Osyka HAT 300 Smart               | 330 090      |
| PK-109        | 2,5 m  | Stockert EP-Shuttle               | 330 079      |
| PK-112        | 2,0 m  | Medtronic Atrak Radionics RFG 3-E | 330 081      |

Distributor:  
**BIOTRONIK GmbH & Co. KG**  
 Woermannkehre 1  
 12359 Berlin - Germany  
 Tel +49 (0)30 68905-0  
 Fax +49 (0)30 6852804  
 sales@biotronik.com



© BIOTRONIK GmbH & Co. KG  
 All rights reserved. Specifications  
 are subject to modification, revision  
 and improvement.

