


CL-ELITE™

VALVED LONG SHEATH INTRODUCER KIT- FEMORAL



- Zero Tolerance Transition Dilator To Sheath, Greatly Reduces Puncture Resistance Locking Dilator To Sheath.
- Smooth Transition From Dilator To Sheath, Greatly Reduces Puncture Resistance Providing Easy And Non-Traumatic Insertion.
- Bends And Flexes Without Gaping or Kinking.
- Minimizes The Possibility Of Damaging Tissue Upon Sheath Entry, Advancement Or Removal. May Lead To Smoother, Easier Vascular Access Even Though Difficult Access Sites.
- Special Designed Lock.
- Silicone Hemostatic Valve.

CONTENTS

- 1- Sheath-Sideport
- 1- Dilator-Locking
- 1- Guidewire
- 1- Angiographic Needle

Ordering Information

CATALOG NO.	DIAMETER	LENGTH (CM)	GUIDEWIRE DIAMETER (INCH)	GUIDEWIRE LENGTH (CM)	NEEDLE GAUGE (GA)
09421	6 F	16cm	0.038"	80	18ga
09422	6 F	23cm	0.038"	80	18ga
09423	7 F	16cm	0.038"	80	18ga
09424	7 F	23cm	0.038"	80	18ga
09425	8 F	16cm	0.038"	80	18ga
09426	8 F	23cm	0.038"	80	18ga

