

CL-ELITE™ VALVED SHEATH INTRODUCER KIT-FEMORAL



- Zero Tolerance Transition Dilator To Sheath, Greatly Reduces Puncture Resistance Locking Dilator To Sheath.
- Smooth Transition From Dilator To Sheath, Greatly Reduces Puncture Resistance Providing Easy And Non-Traumatic Insertion.
- Bends And Flexes Without Gaping or Kinking.
- Minimizes The Possibility Of Damaging Tissue Upon Sheath Entry, Advancement Or Removal. May Lead To Smoother, Easier Vascular Access Even Though Difficult Access Sites.
- Special Designed Lock.
- Silicone Hemostatic Valve.

CONTENTS

- 1- Sheath-Sideport
- 1- Dilator-Locking
- 1- Guidewire
- 1- Angiographic Needle

Ordering Information

CATALOG NO.	DIAMETER	LENGTH (CM)	GUIDEWIRE DIAMETER (INCH)	GUIDEWIRE LENGTH (CM)	NEEDLE GAUGE (GA)
09989	4 F	11cm	0.035"	45	18ga
09990	5 F	11cm	0.035"	45	18ga
09310	5 F	11cm	0.038"	45	18ga
09991	6 F	11cm	0.035"	45	18ga
09312	6 F	11cm	0.038"	45	18ga
09992	7 F	11cm	0.035"	45	18ga
09313	7 F	11cm	0.038"	45	18ga
09993	8 F	11cm	0.035"	45	18ga
09314	8 F	11cm	0.038"	45	18ga
09315	9 F	11cm	0.038"	45	18ga