

ClearView™ MANIFOLD



- RIGHT ROTATING MALE LUER PROVIDES CONVENIENT AND FAST CONNECTIONS
- 360° ROTATING HANDLES AND LUBRICATION ENHANCE HANDLE ROTATION
- VISIBLE HANDLE MAKES THE FLOW CONTROL EASY AND ACCURATE
- PORTS: ONE TO FOUR
- PRESSURE: 250psi or 500 psi
- EXTENSION TUBING: WITH OR WITHOUT
- ON VALVE AND OFF VALVE AVAILABLE

Ordering Information

CATALOG NO.	NO. OF PORTS	DESCRIPTION
14200	2	500 PSI, BLOCK, RIGHT ON HANDLE
14208	2	500 PSI, BLOCK, RIGHT OFF HANDLE
14225	3	500 PSI, BLOCK, RIGHT ON HANDLE
14233	3	500 PSI, BLOCK, RIGHT OFF HANDLE
14229	3	500 PSI, HALF BLOCK, RIGHT ON HANDLE
14237	3	500 PSI, HALF BLOCK, RIGHT OFF HANDLE