


FALCON™ GUIDING CATHETER



- Large Lumen Guide Balancing Flexibility and Support
- Enhanced Visualization
- Full-Wall Technology Construction
- Full Range of Coronary and Peripheral Curves

Ordering Information

PART NUMBER	DESCRIPTION	FR. SIZE	INNER DIAMETER	OUTER DIAMETER	LENGTH
7833-70	JudkinLeft 3.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7875-10	JudkinLeft 3.5	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7875-20	JudkinLeft 4.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7875-30	JudkinLeft 4.5	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7833-G0	JudkinRight 3.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7874-10	JudkinRight 3.5	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7874-20	JudkinRight 4.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7874-30	JudkinRight 4.5	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7874-40	JudkinRight 5.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7877-10	AmplatzLeft 1.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7877-20	AmplatzLeft 2.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7876-10	AmplatzRight 1.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7876-20	AmplatzRight 2.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7833-80	XBU 3.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7833-H0	XBU 3.5	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7838-50	XBU 3.75	6F	0.071"/1.80mm	0.082"/2.08mm	100cm
7838-E0	XBU 4.0	6F	0.071"/1.80mm	0.082"/2.08mm	100cm