


# Fortress

Reinforced Introducer Sheath System / 4F, 5F, 6F  
Indicated to provide access to lower limb arteries



- Coil reinforced sheath for kink and deformation resistance
- Radiopaque sheath marker to provide optimal radiopacity
- Hydrophobic coated sheath<sup>1</sup> for easy insertion
- Smooth tapered sheath and dilator tip

<sup>1</sup> 5F and 6F sheath

# Fortress

## Reinforced Introducer Sheath System / 4F, 5F, 6F

Technical Data		Introducer Sheath	
French size		4F, 5F and 6F	
Length		45 and 100 cm	
Sheath material		Stainless steel coil reinforced polymer, PTFE liner	
Sheath tip		Smooth taper with radiopaque marker band	
Sheath shape		Pre-curved 45 cm Straight 45 and 100 cm	
Sheath coating		Hydrophobic, distal 30 cm (5F and 6F only)	
		Dilator	
Recommended guide wire		0.035"	
Tip		Smooth taper with radiopaque filler material	
Shape		Pre-curved and straight	
Hub		4F: fixed hemostatic valve 5F, 6F: removable hemostatic valve Side-arm with color coded 3-way stop cock	

Ordering Information		Sheath length (cm)	Sheath shape	Introducer description	Catalogue number	Dilator(s) included	Units of issue
4F	45	Curved	Crossover	358813	Curved/Straight	5	
	45	Straight	Straight-45	358814	Curved/Straight	5	
	100	Straight	Straight-100	358815	Curved / Straight	5	
5F	45	Curved	Crossover	386590	Straight	5	
	45	Straight	Straight-45	386591	Straight	5	
	100	Straight	Straight-100	386592	Curved/Straight	5	
6F	45	Curved	Crossover	386593	Straight	5	
	45	Straight	Straight-45	386594	Straight	5	
	100	Straight	Straight-100	386595	Curved/Straight	5	

Fortress is part of the BIOTRONIK **4F** and **6F** Solutions portfolio, including:

- 4F** ■ Guide Wires: **XT-14, Cruiser, Cruiser-18** ■ Balloons: **Passeo-14, Passeo-18**
- Stents: **PRO-Kinetic Energy Explorer, Pulsar-18**
- 6F** ■ Balloons: **Passeo-35, Passeo-35 HP, AngioSculpt** ■ Stents: **Dynamic, Astron, Pulsar-35**

For ordering please contact your  
local sales representative

Manufacturer:  
Contract Medical  
International GmbH  
Lauensteiner Strasse 37  
01277 Dresden · Germany

Distributor:  
BIOTRONIK AG  
Ackerstrasse 6  
8180 Bülach · Switzerland  
Tel +41 (0) 44 8645111  
Fax +41 (0) 44 8645005  
info.vi@biotronik.com  
www.biotronik.com