

EP

Electrophysiology

Diagnostic Catheters


MultiCath

MultiCath

Multipolar catheters for
electrophysiological diagnostics




- Always the right curve shape and pole configuration for each application
- 6F, 5F and 4F versions
- Also available with atraumatic SoftTip



MultiCath 6F

- Very wide range of application-specific curves
- Also available with particularly atraumatic SoftTip


MultiCath 6F	MultiCath 2J 5 mm	MultiCath 2D 5 mm	MultiCath 4J SoftTip 5 mm	MultiCath 4D SoftTip 5 mm
Number of electrodes	2	2	4	4
Electrode spacing	5 mm	5 mm	5 mm	5 mm
Curve shape	Josephson	Damato	Josephson	Damato
Catheter length	110 cm	110 cm	110 cm	110 cm
Patient cable	MPK-4	MPK-4	MPK-4	MPK-4
Order number	351 187	351 190	351185 (MC 4) 351182 (MC 4 SoftTip)	351189 (MC 4) 351181 (MC 4 SoftTip)



MultiCath 5F

- Perfectly suitable for all electrophysiological examinations
- Especially flexible catheter shaft for good tissue contact


MultiCath 5F	MultiCath 4J sp 2 mm	MultiCath 4J sp 5 mm	MultiCath 4J sp 2-5-2 mm	MultiCath 4D 2-5-2 mm
Number of electrodes	4	4	4	4
Electrode spacing	2 mm	5 mm	2-5-2 mm	2-5-2 mm
Curve shape	Josephson special	Josephson special	Josephson special	Damato
Catheter length	80 cm	100 cm	80 cm / 100 cm	80 cm
Patient cable	MPK-4	MPK-4	MPK-4	MPK-4
Order number	351 170	356 434	351 171 (80 cm) 356 433 (100 cm)	351 175




MultiCath 4F

- Enables simple access via brachial vein or jugular vein
- Excellent for pediatric applications


MultiCath 4F	MultiCath 4J sp 5 mm	MultiCath 4J sp 2-5-2 mm	MultiCath 4C 5 mm	MultiCath 4C 2-5-2 mm
Number of electrodes	4	4	4	4
Electrode spacing	5 mm	2-5-2 mm	5 mm	2-5-2 mm
Curve shape	Josephson special	Josephson special	Cournand	Cournand
Catheter length	110 cm	110 cm	110 cm	110 cm
Patient cable	MPK-4	MPK-4	MPK-4	MPK-4
Order number	356 728	356 727	356 726	356 725



MultiCath 4C SoftTip 5 mm	MultiCath 4M SoftTip 5 mm	MultiCath 10J 2 mm	MultiCath 10J 5 mm	MultiCath 10 CS 2-5-2 mm	MultiCath 10 CS 2-8-2 mm
4	4	10	10	10	10
5 mm	5 mm	2 mm	5 mm	2-5-2 mm	2-8-2 mm
Cournand	Multipurpose	Josephson	Josephson	Josephson	Josephson
110 cm	110 cm	110 cm	110 cm	110 cm	110 cm
MPK-4	MPK-4	MPK-10	MPK-10	MPK-10	MPK-10
351 188 (MC 4)	351 194 (MC 4)	351 186	351 193	351 192	351 191
351 183 (MC 4 SoftTip)	351 180 (MC 4 SoftTip)				



MultiCath 4C 2-5-2 mm	MultiCath 4HIS 5 mm	MultiCath 4M 5 mm	MultiCath 6J sp 2-5-2 mm	MultiCath 10J sp 2-5-2 mm	MultiCath 10J sp 2-8-2 mm
4	4	4	6	10	10
2-5-2 mm	5 mm	5 mm	2-5-2 mm	2-5-2 mm	2-8-2 mm
Cournand	HIS	Multipurpose	Josephson special	Josephson special	Josephson special
80 cm	80 cm	80 cm	80 cm	80 cm	80 cm
MPK-4	MPK-4	MPK-4	MPK-10	MPK-10	MPK-10
351 177	359 385	359 387	351 173	359 388	359 389



MultiCath Multipolar diagnostic catheters

Various curve shapes and different electrode spacings provide great flexibility upon selecting a suitable MultiCath catheter for any required application.

- Excellent torqueability
- Outstanding signal quality
- The right curve shape and pole configuration for each application
- Fast and stable positioning

Order Information and Technical Data MultiCath

MultiCath 6F	2-pole		4-pole			
	5 mm	5 mm	5 mm	5 mm	5 mm	5 mm
	Josephson	Damato	Josephson	Damato	Cournand	Multipurpose
	351 187	351 190	351 185	351 189	351 188	351 194
			351 182 (Soft)	351 181 (Soft)	351 183 (Soft)	351 180 (Soft)

10-pole						
	2 mm	5 mm	2-5-2 mm	2-8-2 mm		
	Josephson	Josephson	Josephson	Josephson		
	351 186	351 193	351 192	351 191		

MultiCath 5F	4-pole						
	2 mm	5 mm	2-5-2 mm	2-5-2 mm	2-5-2 mm	5 mm	5 mm
	Josephson sp	Josephson sp	Josephson sp	Damato	Cournand	HIS	Multipurpose
	351 170	356 434	351 171 (80 cm)	351 175	351 177	359 385	359 387
			356 433 (100 cm)				

6-pole		10-pole	
	2-5-2 mm	2-5-2 mm	2-8-2 mm
	Josephson sp	Josephson sp	Josephson sp
	351 173	359 388	359 389


MultiCath 4F	4-pole			
	5 mm	2-5-2 mm	5 mm	2-5-2 mm
	Josephson sp	Josephson sp	Cournand	Cournand
	356 728	356 727	356 726	356 725

Patient cables	MPK-4	MPK-4-R*	MPK-10	MPK-10-R*
Connector device sided	4 × 2 mm pins (touch proof)	4 × 2 mm pins (touch proof)	10 × 2 mm pins (touch proof)	10 × 2 mm pins (touch proof)
Length	2,0 m	2,0 m	2,0 m	2,0 m
Compatible catheters	all 2- and 4-pole MultiCath	all 2- and 4-pole MultiCath	all 6- and 10-pole MultiCath	all 6- and 10-pole MultiCath
Order number	335 805	353 177	335 804	353 173

*-R: autoclavable up to 134° Celsius

Distributor:
 BIOTRONIK GmbH & Co. KG
 Woermannkehre 1
 12359 Berlin · Germany
 Tel +49 (0) 30 68905-0
 Fax +49 (0) 30 6852804
 info@biotronik.com
 www.biotronik.com

999115/H/905



© MultiCath is a registered trademark of
 BIOTRONIK GmbH & Co. KG, Berlin, Germany.

