

**Perfecting** Performance



**Occlutech® Sizing Balloon** Changing the Dimension



# Occlutech® Sizing Balloon

## Changing the Dimension


The Occlutech Sizing Balloon is intended for venous access. It is used without a guide catheter and is introduced directly via the guide wire.



X-ray markers:  
5 mm measured from left edge to left edge

### Designations:

- (1) Usable length
- (2) X-ray markers
- (3) Soft tip
- (4) Inflation port



## Ordering information (dimensions and sizes)

### Occlutech® Sizing Balloon

Article No.	Product Name	Balloon Ø Nominal (mm)	Nominal Fill Vol. (ml)	Max. Fill Vol. (ml)	Balloon Ø Max. Fill Vol. (mm)	Balloon Length (mm)	Usable Length (mm)	Guide Wire (inch)
53SB254	Sizing Balloon	25	25	35	27	45	700	0.035
53SB355	Sizing Balloon	35	60	90	40	50	700	0.035

### Customer Service International

Occlutech, La Cours Gata 2, S-252 31 Helsingborg, Sweden  
Tel +46 (0)42 33 65 32, Fax +46 (0)42 311 09 70  
order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)