

PointSet™ ANGIOGRAPHIC NEEDLE



- EXTREMELY LOW INSERTION FORCES
- REDUCED TISSUE TRAUMA
- TAPERED PINK OR WHITE HUB
- INNER DIAMETER FACILITATES INSERTION OF THE GUIDEWIRE

Ordering Information

| CATALOG NO. | DESCRIPTION | ACCEPTS GW |
|-------------|--|---------------|
| 01316 | PointSET™ Angiographic Needle, 18 Ga. x Thin Wall, 7.0cm, Large Hub, Short Tip, Pink Hub | 0.032"-0.038" |
| 01320 | PointSET™ Angiographic Needle, 18G, 7.0cm, Ultra-Thin Wall, Large Hub With Wing, Long Tip | 0.032"-0.038" |
| 01325 | PointSET™ Angiographic Needle 18G., 7.0cm, Ultra-Thin Wall, Large Hub, Long Tip, Ultra Sharp | 0.032"-0.038" |
| 01330 | PointSET™ Angiographic Needle 21 Ga. x 7cm | 0.032"-0.038" |
| 01352 | PointSET™ Angiographic Needle, 18 Ga. x 2.75" 7.0cm, White Hub | 0.032"-0.038" |