

RadialSTAT-HC™

Radial Sheath Introducer Kits With Hydrophilic Coating



- Zero Tolerance Transition Dilator To Sheath, Greatly Reduces Puncture Resistance Locking Dilator To Sheath.
- Smooth Transition From Dilator To Sheath, Greatly Reduces Puncture Resistance Providing Easy And Non-Traumatic Insertion.
- Bends And Flexes Without Gaping or Kinking.
- Minimizes The Possibility Of Damaging Tissue Upon Sheath Entry, Advancement Or Removal. May Lead To Smoother, Easier Vascular Access Even Though Difficult Access Sites.
- Special Designed Lock. Silicone Hemostatic Valve.

Ordering Information

Part No.	French Size	Length(cm)	Needle Gauge(Ga)	Guidewire Diameter(inch)	Guidewire Length(cm)
06660HC	4F	7	20	0.018	45
06661HC	4F	7	20	0.021	45
06662HC	5F	7	20	0.018	45
06663HC	5F	7	20	0.021	45
06664HC	5F	7	20	0.025	45
06665HC	6F	7	20	0.018	45
06666HC	6F	7	20	0.021	45
06667HC	6F	7	20	0.025	45
06668HC	4F	11	20	0.021	45
06669HC	4F	24	20	0.021	70
09210HC	5F	11	20	0.021	45
09215HC	5F	24	20	0.021	70
09212HC	6F	11	20	0.021	45
09216HC	6F	24	20	0.021	70
09218HC	5F	11	20	0.018	45
09219HC	6F	11	20	0.018	45