

Y-View™ S

PTCA SCREWTOP Y CONNECTOR KITS



- FOR USE ON THE CATHETER PROXIMAL END CONNECTION TO SECURE AND STABILIZE CATHETER PLACEMENT
- MEDICAL GRADE POLYCARBONATE WITH CLOSE TOLERANCE FITTINGS
- 500 psi RATED

Ordering Information

CATALOG NO.	DESCRIPTION
30220	Single Y Connector - Screw Top
30221	Single Y Connector Kit: 1 - Single Y Connector - Screw Top 1-Torque Device 1-Insertion Tool
30233	Single Y Extension Kit: 1-Single Y Connector - Screw Top 1-20cm Extension Line with attached Stopcock 1-Torque Device 1-Insertion Tool

Y-View™ P

PTCA PUSH/PULL OPERATION Y CONNECTOR KIT



- FOR USE ON THE CATHETER PROXIMAL END CONNECTION TO SECURE AND STABILIZE CATHETER PLACEMENT
- MEDICAL GRADE POLYCARBONATE WITH CLOSE TOLERANCE FITTINGS
- 500 psi RATED

Ordering Information

CATALOG NO.	DESCRIPTION
30225	Single Y Connector - Push/Pull
32221	Single Y Connector Kit : 1-Single Y Connector - Push Pull 1-Torque Device 1-Insertion Tool
32640	Single Y Extension Kit : 1-Single Y Connector - Push Pull 1-20cm Extension Line with attached HP Stopcock 1-Torque Device 1-Insertion Tool