

Non-wired Arterial Cannula

Thin Wall

Transparent and Flexible Body

Accurate Positioning in the Artery



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Non-wired Arterial Cannula is made of transparent polyvinyl chloride (PVC), comprising a tip, body and distal end. The tip is provided in straight and curved shapes; it may have a flange for aiding the attachment. The body has a thin and flexible wall with marks for an accurate positioning in the artery. The distal end may have a connector with a luer lock adapter for removal of samples and measurement of blood pressure. Each package contains a sterile and pyrogen-free unit.

Indication: The Wired Arterial Cannula is intended for cannulation of arteries to the bypass cardiopulmonary perfusion system, providing the return of blood to the patient.



| Description | Diameter | Length | Connection | Connector | Code |
|---|----------|--------|------------|----------------------|--------|
| Non-wired Arterial Cannula (Straight Tip) | 08 FR | 16 FR | 1/4" | - | 235778 |
| Non-wired Arterial Cannula (Straight Tip) | 10 FR | 16 FR | 1/4" | - | 235787 |
| Non-wired Arterial Cannula (Straight Tip) | 12 FR | 16 FR | 1/4" | - | 235796 |
| Non-wired Arterial Cannula (Straight Tip) | 14 FR | 16 FR | 1/4" | - | 235802 |
| Non-wired Arterial Cannula (Straight Tip) | 16 FR | 16 FR | 1/4" | - | 235811 |
| Non-wired Arterial Cannula (Straight Tip) | 16 FR | 16 FR | 3/8" | - | 235820 |
| Non-wired Arterial Cannula (Straight Tip) | 18 FR | 16 FR | 3/8" | - | 235839 |
| Non-wired Arterial Cannula (Straight Tip) | 20 FR | 16 FR | 3/8" | - | 235848 |
| Non-wired Arterial Cannula (Straight Tip) | 22 FR | 16 FR | 3/8" | - | 235857 |
| Non-wired Arterial Cannula (Straight Tip) | 24 FR | 16 FR | 3/8" | - | 235866 |
| Non-wired Arterial Cannula (Curved Tip) With Flange | 20 FR | 18 FR | 3/8" | 3/8" x 3/8" com luer | 210720 |
| Non-wired Arterial Cannula (Curved Tip) With Flange | 22 FR | 18 FR | 3/8" | 3/8" x 3/8" com luer | 210739 |
| Non-wired Arterial Cannula (Curved Tip) With Flange | 24 FR | 18 FR | 3/8" | 3/8" x 3/8" com luer | 210748 |
| Non-wired Arterial Cannula (Curved Tip) With Flange | 20 FR | 18 FR | 3/8" | 3/8" x 3/8" com luer | 132851 |
| Non-wired Arterial Cannula (Curved Tip) With Flange | 22 FR | 18 FR | 3/8" | 3/8" x 3/8" com luer | 132842 |
| Non-wired Arterial Cannula (Curved Tip) With Flange | 24 FR | 18 FR | 3/8" | 3/8" x 3/8" com luer | 131432 |



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I

Zip Code 15091-450 • São José do Rio Preto • SP • Brazil

Phone #55 17 2136 7000 • Fax #55 17 2136 7040

Federal Tax ID: 52.828.936/0001-09

Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010

ANVISA Reg. nº 10159030071

www.braile.com.br



Wired Arterial Cannula

Thin Wall

Transparent and Flexible Body

Accurate Positioning in the Artery

Kink Resistant



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Wired Arterial Cannula is made of transparent polyvinyl chloride (PVC), comprising a tip, body and distal end. The tip is provided in straight and curved shapes; it may have a flange for aiding the attachment. The body has a flexible wall with a stainless steel wire reinforcement which avoids kinks while handling it, besides the mark for an accurate positioning in the artery. The distal end may have a connector with a luer lock adapter for removal of samples and measurement of blood pressure. Each package contains a sterile and pyrogen-free unit.

Indication: The Wired Arterial Cannula is intended for cannulation of arteries to the bypass cardiopulmonary perfusion system, providing the patient return of blood.

| Description | Diameter | Length | Connection | Connector | Code |
|--|----------|--------|------------|----------------------|--------|
| Wired Arterial Cannula (Straight Tip) With Flange | 18 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 197018 |
| Wired Arterial Cannula (Straight Tip) With Flange | 20 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 197027 |
| Wired Arterial Cannula (Straight Tip) With Flange | 22 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 197045 |
| Wired Arterial Cannula (Straight Tip) With Flange | 24 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 197054 |
| Wired Arterial Cannula (Curved Tip) 45º | 18 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 196972 |
| Wired Arterial Cannula (Curved Tip) 45º | 20 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 196981 |
| Wired Arterial Cannula (Curved Tip) 45º | 22 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 196990 |
| Wired Arterial Cannula (Curved Tip) 45º | 24 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 197009 |
| Wired Arterial Cannula (Curved Tip) 80º | 18 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 192509 |
| Wired Arterial Cannula (Curved Tip) 80º | 20 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 192387 |
| Wired Arterial Cannula (Curved Tip) 80º | 22 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 192396 |
| Wired Arterial Cannula (Curved Tip) 80º | 24 FR | 30cm | 3/8" | 3/8" x 3/8" com luer | 192402 |
| Wired Arterial Cannula (Straight Tip) With Flange | 08 FR | 13cm | 1/4" | - | 155876 |
| Wired Arterial Cannula (Straight Tip) With Flange | 10 FR | 13cm | 1/4" | - | 155885 |
| Wired Arterial Cannula (Straight Tip) With Flange | 12 FR | 13cm | 1/4" | - | 155894 |
| Wired Arterial Cannula (Straight Tip) With Flange | 14 FR | 15cm | 1/4" | - | 155900 |
| Wired Arterial Cannula (Straight Tip) With Flange | 16 FR | 15cm | 1/4" | - | 155919 |
| Wired Arterial Cannula (Straight Tip) With Flange | 18 FR | 15cm | 1/4" | - | 155928 |
| Wired Arterial Cannula (Straight Tip) With Flange | 20 FR | 15cm | 1/4" | - | 155937 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 08 FR | 13cm | 1/4" | - | 157225 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 10 FR | 13cm | 1/4" | - | 157234 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 12 FR | 13cm | 1/4" | - | 157243 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 14 FR | 15cm | 1/4" | - | 157261 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 16 FR | 15cm | 1/4" | - | 157270 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 18 FR | 15cm | 1/4" | - | 157289 |
| Wired Arterial Cannula (Straight Tip) Without Flange | 20 FR | 15cm | 1/4" | - | 157298 |

| Description | Diameter | Length | Connection | Length | Connector | Code |
|---|----------|--------|------------|--------|----------------------|--------|
| Wired Arterial Cannula (Elongated Straight Tip) | 06 FR | 23cm | 1/4" | 28cm | 1/4" x 1/4" com luer | 238874 |
| Wired Arterial Cannula (Elongated Straight Tip) | 08 FR | 25cm | 1/4" | 30cm | 1/4" x 1/4" com luer | 232283 |
| Wired Arterial Cannula (Elongated Straight Tip) | 10 FR | 25cm | 1/4" | 30cm | 1/4" x 1/4" com luer | 232292 |
| Wired Arterial Cannula (Elongated Straight Tip) | 12 FR | 25cm | 1/4" | 30cm | 1/4" x 1/4" com luer | 232308 |
| Wired Arterial Cannula (Elongated Straight Tip) | 14 FR | 25cm | 1/4" | 30cm | 1/4" x 1/4" com luer | 232317 |
| Wired Arterial Cannula (Elongated Straight Tip) | 16 FR | 25cm | 1/4" | 30cm | 1/4" x 1/4" com luer | 232326 |



Braile Biomédica Indústria, Comércio e Representações Ltda.
Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
CEP 15091-450 • São José do Rio Preto • SP • Brasil
Telefone 17 2136 7000 • Fax 17 2136 7040
Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040
CNPJ: 52.828.936/0001-09
Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010
ANVISA Registration no. 10159030070
www.braile.com.br

