

Veloz Guide Wire


Teflon

Stainless steel coated in PTFE.

Extremely flexible and radiopaque distal end. Straight distal end or in J.

Stiffness levels: Standard, Super Stiff or Extra Stiff.

Diameter of 0.035" (0.89 mm).



ANVISA Reg.: nº 10159030099
www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Braille Description	Braille Code	Model	Body Rigidity	Format of Tip	Coating	Ferrule Length	Wire Diameter	Wire Length
Lunderquist Veloz Guide Wire 0.035"x150cm (J3) Extrastiff Teflon	608370	Lunderquist	Extra Stiff	J3	Teflon	5 cm	0.035"	150 cm
Lunderquist Veloz Guide Wire 0.035"x180cm (J3) Extrastiff Teflon	608371	Lunderquist	Extra Stiff	J3	Teflon	5 cm	0.035"	180 cm
Lunderquist Veloz Guide Wire 0.035"x260cm (J3) Extrastiff Teflon	606495	Lunderquist	Extra Stiff	J3	Teflon	5 cm	0.035"	260 cm
Lunderquist Veloz Guide Wire 0.035"x260cm (J6) Extrastiff Teflon	606586	Lunderquist	Extra Stiff	J6	Teflon	5 cm	0.035"	260 cm
Lunderquist Veloz Guide Wire 0.035"x150cm (J3 7cm) Extrastiff Teflon	608372	Lunderquist	Extra Stiff	J3	Teflon	7 cm	0.035"	150 cm
Lunderquist Veloz Guide Wire 0.035"x180cm (J3 7cm) Extrastiff Teflon	608373	Lunderquist	Extra Stiff	J3	Teflon	7 cm	0.035"	180 cm
Lunderquist Veloz Guide Wire 0.035"x260cm (J3 7cm) Extrastiff Teflon	608374	Lunderquist	Extra Stiff	J3	Teflon	7 cm	0.035"	260 cm
Amplatz Veloz Guide Wire 0.035"x150cm (J3) Superstiff Teflon	606580	Seldinger	Super Stiff	J3	Teflon	4,5 cm	0.035"	150 cm
Amplatz Veloz Guide Wire 0.035"x180cm (J3) Superstiff Teflon	606583	Seldinger	Super Stiff	J3	Teflon	4,5 cm	0.035"	180 cm
Amplatz Veloz Guide Wire 0.035"x260cm (J3) Superstiff Teflon	606585	Seldinger	Super Stiff	J3	Teflon	4,5 cm	0.035"	260 cm
Amplatz Veloz Guide Wire 0.038"x260cm (J3) Superstiff Teflon	608384	Seldinger	Super Stiff	Straight	Teflon	4,5 cm	0.038"	260 cm
Amplatz Veloz Guide Wire 0.035"x150cm (Straight) Superstiff Teflon	608381	Seldinger	Super Stiff	Straight	Teflon	4,5 cm	0.035"	150 cm
Amplatz Veloz Guide Wire 0.035"x180cm (Straight) Superstiff Teflon	608382	Seldinger	Super Stiff	Straight	Teflon	4,5 cm	0.035"	180 cm
Amplatz Veloz Guide Wire 0.035"x260cm (Straight) Superstiff Teflon	608383	Seldinger	Super Stiff	Straight	Teflon	4,5 cm	0.035"	260 cm
Amplatz Veloz Guide Wire 0.038"x260cm (Straight) Superstiff Teflon	608380	Seldinger	Super Stiff	Straight	Teflon	4,5 cm	0.038"	260 cm
Seldinger Veloz Guide Wire 0.035"x260cm (J3) Extrastiff Teflon	608385	Seldinger	Extra Stiff	J3	Teflon	4,5 cm	0.035"	260 cm
Seldinger Veloz Guide Wire 0.038"x260cm (J3) Extrastiff Teflon	608386	Seldinger	Extra Stiff	J3	Teflon	4,5 cm	0.032"	260 cm
Seldinger Veloz Guide Wire 0.035"x260cm (Straight) Extrastiff Teflon	608387	Seldinger	Extra Stiff	Straight	Teflon	4,5 cm	0.035"	260 cm
Seldinger Veloz Guide Wire 0.038"x260cm (Straight) Extrastiff Teflon	608388	Seldinger	Extra Stiff	Straight	Teflon	4,5 cm	0.038"	260 cm
Seldinger Veloz Guide Wire 0.032"x150cm (J3) Standard Teflon	608389	Seldinger	Standard	J3	Teflon	4,5 cm	0.032"	150 cm
Seldinger Veloz Guide Wire 0.035"x150cm (J3) Standard Teflon	606499	Seldinger	Standard	J3	Teflon	4,5 cm	0.035"	150 cm
Seldinger Veloz Guide Wire 0.038"x150cm (J3) Standard Teflon	608390	Seldinger	Standard	J3	Teflon	4,5 cm	0.038"	150 cm
Seldinger Veloz Guide Wire 0.032"x180cm (J3) Standard Teflon	608391	Seldinger	Standard	J3	Teflon	4,5 cm	0.032"	180 cm
Seldinger Veloz Guide Wire 0.032"x180cm (J3) Standard Teflon	606582	Seldinger	Standard	J3	Teflon	4,5 cm	0.035"	180 cm
Seldinger Veloz Guide Wire 0.038"x180cm (J3) Standard Teflon	608392	Seldinger	Standard	J3	Teflon	4,5 cm	0.038"	180 cm
Seldinger Veloz Guide Wire 0.032"x260cm (J3) Standard Teflon	608393	Seldinger	Standard	J3	Teflon	4,5 cm	0.032"	260 cm
Seldinger Veloz Guide Wire 0.035"x260cm (J3) Standard Teflon	606500	Seldinger	Standard	J3	Teflon	4,5 cm	0.035"	260 cm
Seldinger Veloz Guide Wire 0.038"x260cm (J3) Standard Teflon	608394	Seldinger	Standard	J3	Teflon	4,5 cm	0.038"	260 cm
Seldinger Veloz Guide Wire 0.032"x150cm (Straight) Standard Teflon	608395	Seldinger	Standard	Straight	Teflon	4,5 cm	0.032"	150 cm
Seldinger Veloz Guide Wire 0.035"x150cm (Straight) Standard Teflon	606579	Seldinger	Standard	Straight	Teflon	4,5 cm	0.035"	150 cm
Seldinger Veloz Guide Wire 0.038"x150cm (Straight) Standard Teflon	608396	Seldinger	Standard	Straight	Teflon	4,5 cm	0.038"	150 cm
Seldinger Veloz Guide Wire 0.032"x180cm (Straight) Standard Teflon	608397	Seldinger	Standard	Straight	Teflon	4,5 cm	0.032"	180 cm
Seldinger Veloz Guide Wire 0.035"x180cm (Straight) Standard Teflon	606581	Seldinger	Standard	Straight	Teflon	4,5 cm	0.035"	180 cm
Seldinger Veloz Guide Wire 0.038"x180cm (Straight) Standard Teflon	608398	Seldinger	Standard	Straight	Teflon	4,5 cm	0.038"	180 cm
Seldinger Veloz Guide Wire 0.032"x260cm (Straight) Standard Teflon	608399	Seldinger	Standard	Straight	Teflon	4,5 cm	0.032"	260 cm
Seldinger Veloz Guide Wire .035"x260cm (Straight) Standard Teflon	606584	Seldinger	Standard	Straight	Teflon	4,5 cm	0.035"	260 cm
Seldinger Veloz Guide Wire 0.038"x260cm (Straight) Standard Teflon	608400	Seldinger	Standard	Straight	Teflon	4,5 cm	0.038"	260 cm



Braille Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I

Zip Code 15091-450 • São José do Rio Preto • SP • Brazil

Phone #55 17 2136 7000 • Fax #55 17 2136 7040

Federal Tax ID: 52.828.936/0001-09

Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010

ANVISA Reg. nº 10159030099

www.braille.com.br



Veloz Plus Guide Wire

Hydrophilic

Nitinol coated in hydrophilic polyurethane. Excellent radiopaque visualization.
Stiffness level: standard. Facilitated navigation. Safe positioning.
0.035 diameter (0.89 mm).



ANVISA Reg.: nº 10159030100
www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Braile Description	Braile Code	Coating	Material of Core	Format of Tip	Stiffness of Tip	Diameter of Wire	Hardness of Body	Length of the Wire
Veloz Plus Guide Wire Hydrophilic 0.035"x150cm (straight tip) Stiff	608401	Hydrophilic	Nitinol	Straight	Standard	0.035"	Rigid (stiff)	150 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x180cm (straight tip) Stiff	608402	Hydrophilic	Nitinol	Straight	Standard	0.035"	Rigid (stiff)	180 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x260cm (straight tip) Stiff	608403	Hydrophilic	Nitinol	Straight	Standard	0.035"	Rigid (stiff)	260 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x150cm (angular tip) Stiff	608404	Hydrophilic	Nitinol	Angular	Standard	0.035"	Rigid (stiff)	150 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x180cm (angular tip) Stiff	608405	Hydrophilic	Nitinol	Angular	Standard	0.035"	Rígido (stiff)	180 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x260cm (angular tip) Stiff	608406	Hydrophilic	Nitinol	Angular	Standard	0.035"	Rígido (stiff)	260 cm
Veloz Guide Wire Hydrophilic 0.035"x150cm (straight tip) Standard	608407	Hydrophilic	Nitinol	Straight	Standard	0.035"	Standard	150 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x180cm (straight tip) Standard	608408	Hydrophilic	Nitinol	Straight	Standard	0.035"	Standard	180 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x260cm (straight tip) Standard	608409	Hydrophilic	Nitinol	Straight	Standard	0.035"	Standard	260 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x150cm (angular tip)	608410	Hydrophilic	Nitinol	Angular	Standard	0.035"	Standard	150 cm
Standard Veloz Guide Wire Plus Hydrophilic 0.035"x180cm (angular tip)	608411	Hydrophilic	Nitinol	Angular	Standard	0.035"	Standard	180 cm
Standard Veloz Guide Wire Plus Hydrophilic 0.035"x260cm (angular tip)	608412	Hydrophilic	Nitinol	Angular	Standard	0.035"	Standard	260 cm
Standard Veloz Guide Wire Plus Hydrophilic 0.035"x150cm (adaptable tip)	608413	Hydrophilic	Nitinol	Adaptable	Standard	0.035"	Standard	150 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x180cm (adaptable tip)	608414	Hydrophilic	Nitinol	Adaptable	Standard	0.035"	Standard	180 cm
Veloz Plus Guide Wire Hydrophilic 0.035"x260cm (adaptable tip)	608466	Hydrophilic	Nitinol	Adaptable	Standard	0.035"	Standard	260 cm



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 Zip Code 15091-450 • São José do Rio Preto • SP • Brazil
 Phone #55 17 2136 7000 • Fax #55 17 2136 7040
 Federal Tax ID: 52.828.936/0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 ANVISA Reg. nº 10159030100
www.braile.com.br

