

ili[®] INT MEDICAL

ili[®] INT MEDICAL

BEYONDER[®]
Semi Compliant Balloon Catheter

Shanghai Kindly Medical Instruments Co., Ltd.



NO. 925 Jinyuan Yi Road, Shanghai 201803, China



+86 (021) 59140056




+86 (021) 69000013



201803



www.kdl-int.com



NC Non Compliant
Balloon Catheter
NC Beyonder[®]

BEYONDER[®]

Semi Compliant Balloon Catheter




Parameter:

Nominal Pressure (NP): 8 ATM
 Rated Burst Pressure: 16 ATM
 Usable Catheter Length: 140 cm
 Max Guidewire Compatibility: 0.014"


Proximal Shaft Diameter: 1.9 Fr
 Distal Shaft Diameter: 2.4Fr (1.25,1.50mm),
 2.6Fr (2.00~2.50mm),
 2.7Fr (2.75~4.00mm)



Profile size comparison(mm)



Pushing force comparison(N)




*The test was based on 2mm*20mm balloon

Ordering Information

Diameter\ Length	10 mm	15 mm	20 mm	25 mm	30 mm	40 mm
1.25 mm	S012010-4	S012015-4	S012020-4	S012025-4	S012030-4	None
1.50 mm	S015010-4	S015015-4	S015020-4	S015025-4	S015030-4	S015040-4
2.00 mm	S020010-4	S020015-4	S020020-4	S020025-4	S020030-4	S020040-4
2.25 mm	S022010-4	S022015-4	S022020-4	S022025-4	S022030-4	S022040-4
2.50 mm	S025010-4	S025015-4	S025020-4	S025025-4	S025030-4	S025040-4
2.75 mm	S027010-4	S027015-4	S027020-4	S027025-4	S027030-4	S027040-4
3.00 mm	S030010-4	S030015-4	S030020-4	S030025-4	S030030-4	S030040-4
3.25 mm	S032010-4	S032015-4	S032020-4	S032025-4	S032030-4	S032040-4
3.50 mm	S035010-4	S035015-4	S035020-4	S035025-4	S035030-4	S035040-4
4.00 mm	S040010-4	S040015-4	S040020-4	S040025-4	S040030-4	S040040-4
4.50 mm	S045010-4	S045015-4	S045020-4	S045025-4	S045030-4	S045040-4
5.00 mm	S050010-4	S050015-4	S050020-4	S050025-4	S050030-4	S050040-4

Non Compliant Balloon Catheter

NC Beyonder[®]



Parameter:

Nominal Pressure (NP): 12 ATM
 Rated burst pressure: 22 ATM
 Usable Catheter Length: 140 cm
 Max Guidewire Compatibility: 0.014"

Proximal Shaft Diameter: 2.0 Fr
 Distal Shaft Diameter: 2.7Fr



Distal tip comparison



Beyonder Competitor A Competitor B Competitor C Competitor D

Ordering Information

Diameter\ Length	8 mm	10 mm	12 mm	15 mm	18 mm
2.00 mm	NC200008-1	NC200010-1	NC200012-1	NC200015-1	NC200018-1
2.25 mm	NC225008-1	NC225010-1	NC225012-1	NC225015-1	NC225018-1
2.50 mm	NC250008-1	NC250010-1	NC250012-1	NC250015-1	NC250018-1
2.75 mm	NC275008-1	NC275010-1	NC275012-1	NC275015-1	NC275018-1
3.00 mm	NC300008-1	NC300010-1	NC300012-1	NC300015-1	NC300018-1
3.25 mm	NC325008-1	NC325010-1	NC325012-1	NC325015-1	NC325018-1
3.50 mm	NC350008-1	NC350010-1	NC350012-1	NC350015-1	NC350018-1
3.75 mm	NC375008-1	NC375010-1	NC375012-1	NC375015-1	NC375018-1
4.00 mm	NC400008-1	NC400010-1	NC400012-1	NC400015-1	NC400018-1
4.50 mm	NC450008-1	NC450010-1	NC450012-1	NC450015-1	NC450018-1
5.00 mm	NC500008-1	NC500010-1	NC500012-4	NC500015-1	NC500018-1