



# Extension Tube



- Unique braiding technology provides clear visualization of air bubbles inside the tubing
- Length customization is acceptable

Ref.NO.	Model	Production Description	Ref.NO.	Model	Production Description
PMT0101	30cm	Female-Male Luer Lock No Leakage under 500psi Pressure	PMT0501	30cm	Rotating adapter No Leakage under 1200psi Pressure
PMT0102	60cm		PMT0502	60cm	
PMT0103	90cm		PMT0503	90cm	
PMT0104	120cm		PMT0504	120cm	
PMT0105	150cm		PMT0505	150cm	
PMT0106	180cm		PMT0509	180cm	
PMT0301	30cm	Female-Male Luer Lock No Leakage under 1200psi Pressure	PMT0201	30cm	3-way stopcock-Male Luer lock No Leakage under 40psi Pressure
PMT0302	60cm		PMT0202	60cm	
PMT0303	90cm		PMT0203	90cm	
PMT0304	120cm		PMT0204	120cm	
PMT0304A	150cm		Male-Male Luer Lock No Leakage under 500psi Pressure	PMT0101P	30cm
PMT0304B	180cm				

# Extension Tube



- Unique braiding technology provides clear visualization of air bubbles inside the tubing

Ref.No.	Product Description	Model
PMT0801A	30CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	30cm
PMT0801B	50CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	50cm
PMT0801C	60CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	60cm
PMT0801D	90CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	90cm
PMT0801E	120CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	120cm
PMT0801F	150CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	150cm
PMT0801G	180CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	180cm
PMT0801H	200CM*1200PSI, 1.5*4.7MM, PVC, Female-Male Luer-Lock	200cm
PMT0901A	30CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	30cm
PMT0901B	50CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	50cm
PMT0901C	60CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	60cm
PMT0901D	90CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	90cm
PMT0901E	120CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	120cm
PMT0901F	150CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	150cm
PMT0901G	180CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	180cm
PMT0901H	200CM*1200PSI, 1.5*4.7MM, PVC, Rotating adapter	200cm