


Pressure Transducer



- Consistent and accurate pressure monitoring with pressure range: -4~40kPa (-30~300mmHg)
- With built-in flush device to simplify priming the tube and flushing
- Sensitivity: 5 μ V/V/mmHg \pm 2%

Ref.NO.	Model	Monitor cable connector	Optional accessories	Ref.NO.	Model	Monitor cable connector	Optional accessories
DPT0101	DPT0101	BD	Flush Set	DPT0403	DPT0403	Argon	Pressure monitoring tubing (1200mm), Flush Set
DPT0102	DPT0102		Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)	DPT0404	DPT0404		Flush Set, Infusion set (1500mm)
DPT0103	DPT0103		Pressure monitoring tubing (1200mm), Flush Set	DPT0501	DPT0501	ACE	Flush Set
DPT0104	DPT0104		Flush Set, Infusion set (1500mm)	DPT0502	DPT0502		Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)
DPT0201	DPT0201	Philip	Flush Set	DPT0503	DPT0503		Pressure monitoring tubing (1200mm), Flush Set
DPT0202	DPT0202		Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)	DPT0504	DPT0504	Flush Set, Infusion set (1500mm)	
DPT0203	DPT0203		Pressure monitoring tubing (1200mm), Flush Set	DPT0601	DPT0601	Utah	Flush Set
DPT0204	DPT0204		Flush Set, Infusion set (1500mm)	DPT0602	DPT0602		Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)
DPT0301	DPT0301	Abbot	Flush Set	DPT0603	DPT0603		Pressure monitoring tubing (1200mm), Flush Set
DPT0302	DPT0302		Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)	DPT0604	DPT0604	Flush Set, Infusion set (1500mm)	
DPT0303	DPT0303		Pressure monitoring tubing (1200mm), Flush Set	DPT0701	DPT0701	ELS	Flush Set
DPT0304	DPT0304		Flush Set, Infusion set (1500mm)	DPT0702	DPT0702		Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)
DPT0401	DPT0401	Argon	Flush Set	DPT0703	DPT0703		Pressure monitoring tubing (1200mm), Flush Set
DPT0402	DPT0402		Flush Pressure monitoring tubing (1200mm), Flush Set, Infusion set (1500mm)	DPT0704	DPT0704		Flush Set, Infusion set (1500mm)