




LonAccess™

Guiding Sheath




Hydrophilic coating on distal end

Stainless Steel Coil

Radiopaque marker

PTFE Inner Layer

Guiding Sheath



Ordering information

French Size	REF.NO.	Usable Length(cm)	Hydrophilic Coating Length(cm)	Curve Style	Curve Shape	
5Fr	VSS5F01	40	30	Straight		
	VSS5F02	55	40	Straight		
	VSS5F03	70	50	Straight		
	VSS5F04	90	60	Straight		
	VSM5F01	40	30	Curved		
	VSM5F02	55	40	Curved		
	VSM5F03	70	50	Curved		
	VSM5F04	90	60	Curved		
	VSMP5F01	40	30	Multi-Purpose		
	VSMP5F02	55	40	Multi-Purpose		
	VSMP5F03	70	50	Multi-Purpose		
	VSMP5F04	90	60	Multi-Purpose		
	6Fr	VSS6F01	40	30	Straight	
		VSS6F02	55	40	Straight	
		VSS6F03	70	50	Straight	
		VSS6F04	90	60	Straight	
VSM6F01		40	30	Curved		
VSM6F02		55	40	Curved		
VSM6F03		70	50	Curved		
VSM6F04		90	60	Curved		
VSMP6F01		40	30	Multi-Purpose		
VSMP6F02		55	40	Multi-Purpose		
VSMP6F03		70	50	Multi-Purpose		
VSMP6F04		90	60	Multi-Purpose		

Packed 1 per box and includes dilator. All dilators are 0.038" wire compatible.

Guiding Sheath

Ordering information

French Size	REF.NO.	Usable Length(cm)	Hydrophilic Coating Length(cm)	Curve Style	Curve Shape	
7Fr	VSS7F01	40	30	Straight		
	VSS7F02	55	40	Straight		
	VSS7F03	70	50	Straight		
	VSS7F04	90	60	Straight		
	VSM7F01	40	30	Curved		
	VSM7F02	55	40	Curved		
	VSM7F03	70	50	Curved		
	VSM7F04	90	60	Curved		
	VSMP7F01	40	30	Multi-Purpose		
	VSMP7F02	55	40	Multi-Purpose		
	VSMP7F03	70	50	Multi-Purpose		
	VSMP7F04	90	60	Multi-Purpose		
	8Fr	VSS8F01	40	30	Straight	
		VSS8F02	55	40	Straight	
		VSS8F03	70	50	Straight	
		VSS8F04	90	60	Straight	
VSM8F01		40	30	Curved		
VSM8F02		55	40	Curved		
VSM8F03		70	50	Curved		
VSM8F04		90	60	Curved		
VSMP8F01		40	30	Multi-Purpose		
VSMP8F02		55	40	Multi-Purpose		
VSMP8F03		70	50	Multi-Purpose		
VSMP8F04		90	60	Multi-Purpose		

Packed 1 per box and includes dilator. All dilators are 0.038" wire compatible.