




Angiography Catheter

**POINT ⚡ CURVE**<sup>™</sup>



# POINT CURVE™



## **Pebax and Stainless-Steel braiding**

Provides outstanding pushability and crossability.

## **Double braiding with 32 strands**

Designed to achieve 1:1 torque and for excellent support.

## **Soft tip**

Reduces potential for vessel trauma.

## **Large inner lumen and ultra-thin wall**

Ensures high flow rate of contrast medical.

## **Reinforced shaft**

Can resist up to 1200psi pressure and maintain flexibility.

## **Inner Lumen Sizes**

4F 0.042" (1.07mm)

5F 0.046" (1.17mm)


6F 0.054" (1.37mm)

7F 0.064" (1.62mm)

Color hub for easy identification

Maximum guidewire size: 0.038"/0.97mm

Maximum pressure: 1200psi/83bar




# POINT CURVE™

## CARDIOLOGY CATHETERS


Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
<b>JUDKINS LEFT(JL)</b>				
ACJL354F1100	4	0.042"/1.07mm	100cm	JL3.5*4F
ACJL404F1100	4	0.042"/1.07mm	100cm	JL4.0*4F
ACJL504F1100	4	0.042"/1.07mm	100cm	JL5.0*4F
ACJL604F1100	4	0.042"/1.07mm	100cm	JL6.0*4F
ACJL355F1100	5	0.046"/1.17mm	100cm	JL3.5*5F
ACJL355F1125	5	0.046"/1.17mm	125cm	JL3.5*5F
ACJL405F1100	5	0.046"/1.17mm	100cm	JL4.0*5F
ACJL405F1125	5	0.046"/1.17mm	125cm	JL4.0*5F
ACJL405F1130	5	0.046"/1.17mm	130cm	JL4.0*5F
ACJL455F1100	5	0.046"/1.17mm	100cm	JL4.5*5F
ACJL505F1100	5	0.046"/1.17mm	100cm	JL5.0*5F
ACJL505F1125	5	0.046"/1.17mm	125cm	JL5.0*5F
ACJL605F1100	5	0.046"/1.17mm	100cm	JL6.0*5F
ACJL605F1125	5	0.046"/1.17mm	125cm	JL6.0*5F
ACJL356F1100	6	0.054"/1.37mm	100cm	JL3.5*6F
ACJL356F1125	6	0.054"/1.37mm	125cm	JL3.5*6F
ACJL406F1100	6	0.054"/1.37mm	100cm	JL4.0*6F
ACJL406F1125	6	0.054"/1.37mm	125cm	JL4.0*6F
ACJL456F1100	6	0.054"/1.37mm	100cm	JL4.5*6F
ACJL506F1100	6	0.054"/1.37mm	100cm	JL5.0*6F
ACJL506F1125	6	0.054"/1.37mm	125cm	JL5.0*6F
ACJL606F1100	6	0.054"/1.37mm	100cm	JL6.0*6F
ACJL606F1125	6	0.054"/1.37mm	125cm	JL6.0*6F
ACJL357F1100	7	0.064"/1.62mm	100cm	JL3.5*7F
ACJL407F1100	7	0.064"/1.62mm	100cm	JL4.0*7F
ACJL457F1100	7	0.064"/1.62mm	100cm	JL4.5*7F
ACJL507F1100	7	0.064"/1.62mm	100cm	JL5.0*7F
ACJL607F1100	7	0.064"/1.62mm	100cm	JL6.0*7F




**JL3.5**




**JL4.0**



**JL4.5**



**JL5.0**





**JL6.5**

Recommended guide wire is 0.038" (0.97mm)

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
----------------	-------------	------------------------	-----------	-------

### JUDKINS RIGHT(JR)

ACJR354F1100	4	0.042"/1.07mm	100cm	JR3.5*4F
ACJR404F1100	4	0.042"/1.07mm	100cm	JR4.0*4F
ACJR504F1100	4	0.042"/1.07mm	100cm	JR5.0*4F
ACJR604F1100	4	0.042"/1.07mm	100cm	JR6.0*4F
ACJR355F1100	5	0.046"/1.17mm	100cm	JR3.5*5F
ACJR355F1125	5	0.046"/1.17mm	125cm	JR3.5*5F
ACJR405F1100	5	0.046"/1.17mm	100cm	JR4.0*5F
ACJR405F1125	5	0.046"/1.17mm	125cm	JR4.0*5F
ACJR405F1130	5	0.046"/1.17mm	130cm	JR4.0*5F
ACJR505F1100	5	0.046"/1.17mm	100cm	JR5.0*5F
ACJR505F1125	5	0.046"/1.17mm	125cm	JR5.0*5F
ACJR605F1100	5	0.046"/1.17mm	100cm	JR6.0*5F
ACJR605F1125	5	0.046"/1.17mm	125cm	JR6.0*5F
ACJR356F1100	6	0.054"/1.37mm	100cm	JR3.5*6F
ACJR356F1125	6	0.054"/1.37mm	125cm	JR3.5*6F
ACJR406F1100	6	0.054"/1.37mm	100cm	JR4.0*6F
ACJR406F1125	6	0.054"/1.37mm	125cm	JR4.0*6F
ACJR506F1100	6	0.054"/1.37mm	100cm	JR5.0*6F
ACJR506F1125	6	0.054"/1.37mm	125cm	JR5.0*6F
ACJR606F1100	6	0.054"/1.37mm	100cm	JR6.0*6F
ACJR606F1125	6	0.054"/1.37mm	125cm	JR6.0*6F
ACJR357F1100	7	0.064"/1.62mm	100cm	JR3.5*7F
ACJR407F1100	7	0.064"/1.62mm	100cm	JR4.0*7F
ACJR507F1100	7	0.064"/1.62mm	100cm	JR5.0*7F
ACJR607F1100	7	0.064"/1.62mm	100cm	JR6.0*7F



### CORONARY BYPASS RIGHT(CBR)

ACCBR4F1100	4	0.042"/1.07mm	100cm	CBR*4F
ACCBR5F1100	5	0.046"/1.17mm	100cm	CBR*5F
ACCBR6F1100	6	0.054"/1.37mm	100cm	CBR*6F
ACCBR7F1100	7	0.064"/1.62mm	100cm	CBR*7F



### CORONARY BYPASS LEFT(CBL)

ACCBL 4F1100	4	0.042"/1.07mm	100cm	CBL *4F
ACCBL 5F1100	5	0.046"/1.17mm	100cm	CBL *5F
ACCBL 6F1100	6	0.054"/1.37mm	100cm	CBL *6F
ACCBL 7F1100	7	0.064"/1.62mm	100cm	CBL *7F




Recommended guide wire is 0.038" (0.97mm)

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
----------------	-------------	------------------------	-----------	-------




### AMPLATZ LEFT(AL)

ACAL14F1100	4	0.042"/1.07mm	100cm	AL1*4F
ACAL24F1100	4	0.042"/1.07mm	100cm	AL2*4F
ACAL34F1100	4	0.042"/1.07mm	100cm	AL3*4F
ACAL44F1100	4	0.042"/1.07mm	100cm	AL4*4F
ACAL15F1100	5	0.046"/1.17mm	100cm	AL1*5F
ACAL25F1100	5	0.046"/1.17mm	100cm	AL2*5F
ACAL35F1100	5	0.046"/1.17mm	100cm	AL3*5F
ACAL45F1100	5	0.046"/1.17mm	100cm	AL4*5F
ACAL16F1100	6	0.054"/1.37mm	100cm	AL1*6F
ACAL26F1100	6	0.054"/1.37mm	100cm	AL2*6F
ACAL36F1100	6	0.054"/1.37mm	100cm	AL3*6F
ACAL46F1100	6	0.054"/1.37mm	100cm	AL4*6F
ACAL17F1100	7	0.064"/1.62mm	100cm	AL1*7F
ACAL27F1100	7	0.064"/1.62mm	100cm	AL2*7F
ACAL37F1100	7	0.064"/1.62mm	100cm	AL3*7F
ACAL47F1100	7	0.064"/1.62mm	100cm	AL4*7F



### AMPLATZ RIGHT(AR)

ACAR14F1100	4	0.042"/1.07mm	100cm	AR1*4F
ACAR24F1100	4	0.042"/1.07mm	100cm	AR2*4F
ACAR34F1100	4	0.042"/1.07mm	100cm	AR3*4F
ACAR15F1100	5	0.046"/1.17mm	100cm	AR1*5F
ACAR25F1100	5	0.046"/1.17mm	100cm	AR2*5F
ACAR35F1100	5	0.046"/1.17mm	100cm	AR3*5F
ACAR16F1100	6	0.054"/1.37mm	100cm	AR1*6F
ACAR26F1100	6	0.054"/1.37mm	100cm	AR2*6F
ACAR36F1100	6	0.054"/1.37mm	100cm	AR3*6F
ACAR17F1100	7	0.064"/1.62mm	100cm	AR1*7F
ACAR27F1100	7	0.064"/1.62mm	100cm	AR2*7F
ACAR37F1100	7	0.064"/1.62mm	100cm	AR3*7F



### INTERNAL MARMARY(IM)

ACIM4F1100	4	0.042"/1.07mm	100cm	IM*4F
ACIM5F1100	5	0.046"/1.17mm	100cm	IM*5F
ACIM6F1100	6	0.054"/1.37mm	100cm	IM*6F
ACIM7F1100	7	0.064"/1.62mm	100cm	IM*7F



Recommended guide wire is 0.038" (0.97mm)

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model	Side Holes
----------------	-------------	------------------------	-----------	-------	------------

**PIGTAIL(PIG)**


ACPIGS4F1080	4	0.042"/1.07mm	80cm	PIG*4F	6
ACPIGS4F1110	4	0.042"/1.07mm	110cm	PIG*4F	6
ACPIGS4F1125	4	0.042"/1.07mm	125cm	PIG*4F	6
ACPIGS4F1130	4	0.042"/1.07mm	130cm	PIG*4F	6
ACPIGS5F1065	5	0.046"/1.17mm	65cm	PIG*5F	6
ACPIGS5F1110	5	0.046"/1.17mm	110cm	PIG*5F	6
ACPIGS5F1125	5	0.046"/1.17mm	125cm	PIG*5F	6
ACPIGS5F1130	5	0.046"/1.17mm	130cm	PIG*5F	6
ACPIGS6F1080	6	0.054"/1.37mm	80cm	PIG*6F	5
ACPIGS6F1110	6	0.054"/1.37mm	110cm	PIG*6F	5
ACPIGS6F1125	6	0.054"/1.37mm	125cm	PIG*6F	5
ACPIGS7F1080	7	0.064"/1.62mm	80cm	PIG*7F	5
ACPIGS7F1110	7	0.064"/1.62mm	110cm	PIG*7F	5
ACPIGS7F1125	7	0.064"/1.62mm	125cm	PIG*7F	5



**PIG**

**PIGTAIL 145°(PIG)**

ACPIGA5F1110	5	0.046"/1.17mm	110cm	PIG*5F	6
ACPIGA6F1110	6	0.054"/1.37mm	110cm	PIG*6F	6




**PIG 145°**

**PIGTAIL 155°(PIG)**

ACPIGB5F1110	5	0.046"/1.17mm	110cm	PIG*5F	5
ACPIGB6F1110	6	0.054"/1.37mm	110cm	PIG*6F	5

**ULTIMATERADIAL(TIG)**

ACTIG4F1100	4	0.042"/1.07mm	100cm	TIG*4F	1
ACTIG5F1100	5	0.046"/1.17mm	100cm	TIG*5F	1
ACTIG6F1100	6	0.054"/1.37mm	100cm	TIG*6F	1
ACTIG7F1100	7	0.064"/1.62mm	100cm	TIG*7F	1



**TIG**

Recommended guide wire is 0.038" (0.97mm)

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model	Side Holes
<b>MULTIPURPOSE(MPA)</b>					
ACMPA14F1100	4	0.042"/1.07mm	100cm	MultipurposeA1*4F	0
ACMPA24F1100	4	0.042"/1.07mm	100cm	MultipurposeA2*4F	2
ACMPA15F1100	5	0.046"/1.17mm	100cm	MultipurposeA1*5F	0
ACMPA25F1100	5	0.046"/1.17mm	100cm	MultipurposeA2*5F	2
ACMPA16F1100	6	0.054"/1.37mm	100cm	MultipurposeA1*6F	0
ACMPA26F1100	6	0.054"/1.37mm	100cm	MultipurposeA2*6F	2
ACMPA17F1100	7	0.064"/1.62mm	100cm	MultipurposeA1*7F	0
ACMPA27F1100	7	0.064"/1.62mm	100cm	MultipurposeA2*7F	2



**MPA1**



**MPA2**

**MULTIPURPOSE(MPB)**

ACMPB14F1100	4	0.042"/1.07mm	100cm	MultipurposeB1*4F	0
ACMPB24F1100	4	0.042"/1.07mm	100cm	MultipurposeB2*4F	2
ACMPB15F1100	5	0.046"/1.17mm	100cm	MultipurposeB1*5F	0
ACMPB25F1100	5	0.046"/1.17mm	100cm	MultipurposeB2*5F	2
ACMPB16F1100	6	0.054"/1.37mm	100cm	MultipurposeB1*6F	0
ACMPB26F1100	6	0.054"/1.37mm	100cm	MultipurposeB2*6F	2
ACMPB17F1100	7	0.064"/1.62mm	100cm	MultipurposeB1*7F	0
ACMPB27F1100	7	0.064"/1.62mm	100cm	MultipurposeB2*7F	2



**MPB1**






**MPB2**

Recommended guide wire is 0.038" (0.97mm)

## PERIPHERAL CATHETERS

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
<b>HEADHUNTER(H)</b>				
ACH14F0100	4	0.042"/1.07mm	100cm	Headhunter1*4F
ACH34F0100	4	0.042"/1.07mm	100cm	Headhunter3*4F
ACH15F0100	5	0.046"/1.17mm	100cm	Headhunter1*5F
ACH35F0100	5	0.046"/1.17mm	100cm	Headhunter3*5F
ACH16F0100	6	0.054"/1.37mm	100cm	Headhunter1*6F
ACH36F0100	6	0.054"/1.37mm	100cm	Headhunter3*6F
ACH17F0100	7	0.064"/1.62mm	100cm	Headhunter1*7F
ACH37F0100	7	0.064"/1.62mm	100cm	Headhunter3*7F
<b>COBRA(COB)</b>				
ACCOBRA14F1100	4	0.042"/1.07mm	100cm	COBRA1*4F
ACCOBRA24F1100	4	0.042"/1.07mm	100cm	COBRA2*4F
ACCOBRA15F1100	5	0.046"/1.17mm	100cm	COBRA1*5F
ACCOBRA25F1100	5	0.046"/1.17mm	100cm	COBRA2*5F
ACCOBRA16F1100	6	0.054"/1.37mm	100cm	COBRA1*6F
ACCOBRA26F1100	6	0.054"/1.37mm	100cm	COBRA2*6F
ACCOBRA17F1100	7	0.064"/1.62mm	100cm	COBRA1*7F
ACCOBRA27F1100	7	0.064"/1.62mm	100cm	COBRA2*7F
<b>HINCK</b>				
ACHINCK4F1100	4	0.042"/1.07mm	100cm	HINCK*4F
ACHINCK5F1100	5	0.046"/1.17mm	100cm	HINCK*5F
ACHINCK6F1100	6	0.054"/1.37mm	100cm	HINCK*6F
ACHINCK7F1100	7	0.064"/1.62mm	100cm	HINCK*7F



Recommended guide wire is 0.038" (0.97mm)



Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
----------------	-------------	------------------------	-----------	-------

**BENTSON(SON)**

ACSONI4F1100	4	0.042"/1.07mm	100cm	SON I*4F
ACSONIW4F1100	4	0.042"/1.07mm	100cm	SON IW*4F
ACSONII4F1100	4	0.042"/1.07mm	100cm	SON II*4F
ACSONI IW4F1100	4	0.042"/1.07mm	100cm	SON IIW*4F
ACSONIII4F1100	4	0.042"/1.07mm	100cm	SON III*4F
ACSONIIIW4F1100	4	0.042"/1.07mm	100cm	SON IIIW*4F

ACSONI5F1100	5	0.046"/1.17mm	100cm	SON I*5F
ACSONIW5F1100	5	0.046"/1.17mm	100cm	SON IW*5F
ACSONII5F1100	5	0.046"/1.17mm	100cm	SON II*5F
ACSONI IW5F1100	5	0.046"/1.17mm	100cm	SON IIW*5F
ACSONIII5F1100	5	0.046"/1.17mm	100cm	SON III*5F
ACSONIIIW5F1100	5	0.046"/1.17mm	100cm	SON IIIW*5F

ACSONI6F1100	6	0.054"/1.37mm	100cm	SON I*6F
ACSONIW6F1100	6	0.054"/1.37mm	100cm	SON IW*6F
ACSONII6F1100	6	0.054"/1.37mm	100cm	SON II*6F
ACSONI IW6F1100	6	0.054"/1.37mm	100cm	SON IIW*6F
ACSONIII6F1100	6	0.054"/1.37mm	100cm	SON III*6F
ACSONIIIW6F1100	6	0.054"/1.37mm	100cm	SON IIIW*6F

ACSONI7F1100	7	0.064"/1.62mm	100cm	SON I*7F
ACSONIW7F1100	7	0.064"/1.62mm	100cm	SON IW*7F
ACSONII7F1100	7	0.064"/1.62mm	100cm	SON II*7F
ACSONI IW7F1100	7	0.064"/1.62mm	100cm	SON IIW*7F
ACSONIII7F1100	7	0.064"/1.62mm	100cm	SON III*7F
ACSONIIIW7F1100	7	0.064"/1.62mm	100cm	SON IIIW*7F

**MIKAELSSON(MP)**

ACMP4F1080	4	0.042"/1.07mm	80cm	MP*4F
------------	---	---------------	------	-------

ACMP5F1080	5	0.046"/1.17mm	80cm	MP*5F
------------	---	---------------	------	-------

ACMP6F1080	6	0.054"/1.37mm	80cm	MP*6F
------------	---	---------------	------	-------

ACMP7F1080	7	0.064"/1.62mm	80cm	MP*7F
------------	---	---------------	------	-------





Recommended guide wire is 0.038" (0.97mm)

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
----------------	-------------	------------------------	-----------	-------

**SIMMONS(SIM)**

ACSIM14F1100	4	0.042"/1.07mm	100cm	Simmons1*4F
ACSIM24F1100	4	0.042"/1.07mm	100cm	Simmons2*4F
ACSIM34F1100	4	0.042"/1.07mm	100cm	Simmons3*4F
ACSIM15F1100	5	0.046"/1.17mm	100cm	Simmons1*5F
ACSIM25F1100	5	0.046"/1.17mm	100cm	Simmons2*5F
ACSIM25F1125	5	0.046"/1.17mm	125cm	Simmons2*5F
ACSIM35F1100	5	0.046"/1.17mm	100cm	Simmons3*5F
ACSIM35F1125	5	0.046"/1.17mm	125cm	Simmons3*5F
ACSIM16F1100	6	0.054"/1.37mm	100cm	Simmons1*6F
ACSIM26F1100	6	0.054"/1.37mm	100cm	Simmons2*6F
ACSIM36F1100	6	0.054"/1.37mm	100cm	Simmons3*6F
ACSIM17F1100	7	0.064"/1.62mm	100cm	Simmons1*7F
ACSIM27F1100	7	0.064"/1.62mm	100cm	Simmons2*7F
ACSIM37F1100	7	0.064"/1.62mm	100cm	Simmons3*7F



**DAVIS**

ACDAVIS4F1100	4	0.042"/1.07mm	100cm	DAVIS*4F
ACDAVIS5F1100	5	0.046"/1.17mm	100cm	DAVIS*5F
ACDAVIS6F1100	6	0.054"/1.37mm	100cm	DAVIS*6F
ACDAVIS7F1100	7	0.064"/1.62mm	100cm	DAVIS*7F



**RENAL**

ACRENAL4F1100	4	0.042"/1.07mm	100cm	RENAL*4F
ACRENAL5F1100	5	0.046"/1.17mm	100cm	RENAL*5F
ACRENAL6F1100	6	0.054"/1.37mm	100cm	RENAL*6F
ACRENAL7F1100	7	0.064"/1.62mm	100cm	RENAL*7F




Recommended guide wire is 0.038" (0.97mm)

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model
----------------	-------------	------------------------	-----------	-------

**HEPATICDUCTR(RH)**


ACRH4F1100	4	0.042"/1.07mm	100cm	RH*4F
ACRH5F1100	5	0.046"/1.17mm	100cm	RH*5F
ACRH6F1100	6	0.054"/1.37mm	100cm	RH*6F
ACRH7F1100	7	0.064"/1.62mm	100cm	RH*7F



**RH**

**YASHIRO**

ACYASHIRO4F1100	4	0.042"/1.07mm	100cm	Yashiro*4F
ACYASHIRO5F1100	5	0.046"/1.17mm	100cm	Yashiro*5F
ACYASHIRO6F1100	6	0.054"/1.37mm	100cm	Yashiro*6F
ACYASHIRO7F1100	7	0.064"/1.62mm	100cm	Yashiro*7F



**YASHIRO**

Catalog Number	French Size	Inner Diameter inch/mm	Length cm	Model	Side Holes	Radiopaque Marker
----------------	-------------	------------------------	-----------	-------	------------	-------------------

**MULTIPURPOSE(MPA)**

ACPIG14F52M65	4	0.042"/1.07mm	65cm	PIG*4F	5	2*(2cm apart)
ACPIG14F52M100	4	0.042"/1.07mm	100cm	PIG*4F	5	2*(2cm apart)
ACPIG4F520M65	4	0.042"/1.07mm	65cm	PIG*4F	5	20*(1cm apart)
ACPIG4F520M100	4	0.042"/1.07mm	100cm	PIG*4F	5	20*(1cm apart)
ACPIG15F102M65	5	0.046"/1.17mm	65cm	PIG*5F	10	2*(2cm apart)
ACPIG15F102M100	5	0.046"/1.17mm	100cm	PIG*5F	10	2*(2cm apart)
ACPIG5F520M65	5	0.046"/1.17mm	65cm	PIG*5F	5	20*(1cm apart)
ACPIG5F520M100	5	0.046"/1.17mm	100cm	PIG*5F	5	20*(1cm apart)
ACPIG5F1020M65	5	0.046"/1.17mm	65cm	PIG*5F	10	20*(1cm apart)
ACPIG5F1020M100	5	0.046"/1.17mm	100cm	PIG*5F	10	20*(1cm apart)



**PIG**

Recommended guide wire is 0.038" (0.97mm)



**Add**

No.925 Jinyuan Rd.1,  
Shanghai China



**TEL**

021-59140056



**FAX**

021-69000013



**P.C.**

201803



**WEB**

www.kdl-int.com



**Shanghai Kindly Medical Instruments Co., Ltd.**