


P-14TM & P-18TM

DenaliWireTM Guidewire



Shanghai Kindly Medical Instruments Co., Ltd.

P-14™ & P-18™

DenaliWire™ Guidewire

Designed for navigating and crossing stenosed lesions above and below in peripheral.

P-14™ DenaliWire™

REF.NO.	Diameter(Inch)	Length(cm)	Polymer Length(cm)	Tip Shape
HG144SF	0.014	180	8	Straight
HG148SF	0.014	300	8	Straight

P-18™ DenaliWire™

REF.NO.	Diameter(Inch)	Length(cm)	Polymer Length(cm)	Tip Shape
HG184SF	0.018	180	8	Straight
HG188SF	0.018	300	8	Straight



Shanghai Kindly Medical Instruments Co., Ltd.

📍 Add: No.925 Jinyuan Rd 1, Shanghai, China

☎ TEX: +86-21-59140056

📠 FAX: +86-21-69000013

✉ P.C.: 201803