


# GUIDE WIRE

## Peripheral Vascular



### Parameters

Total length: 185cm/300cm

Tip radiopacity: 3.0cm

Tip Shape: Straight / J

Core material: Linear elastic nitinol

### \*FEATURES

#### Dual core technology

For smooth transition between Nitinol core to SS304V core.

#### SS304V core with PTFE coating

Providing improved smooth device delivery and guide wire trackability.


#### Tungsten based polymer jacket with hydrophilic coating

Enables enhanced visualization and steering ability.

#### Distal Nitinol design

For excellent durability and tip shape retention.

# Peripheral Vascular Guide Wire



- Sina(0.014"-0.5gf)
- Sina(0.014"-0.7gf)
- BW (0.014"-0.6gf)
- BW (0.014"-0.8gf)
- BW (0.014"-1.0gf)
- PV3 (0.014"-2.5gf)
- PV3 (0.014"-3.0gf)
- PV3 Plus (0.018"-2.5gf)
- PV3 Plus (0.018"-3.0gf)



## ORDERING INFORMATION

REF.NO.	Diameter	Length	Tip Load	Spring coil length	Radiopaque length	Tip Shape
NS140303BS25	0.014 inch	185cm	2.5gf	3.0cm	3cm	Straight
NS140303BJ25	0.014 inch	185cm	2.5gf	3.0cm	3cm	J
NS140303FS25	0.014 inch	300cm	2.5gf	3.0cm	3cm	Straight
NS140303FJ25	0.014 inch	300cm	2.5gf	3.0cm	3cm	J
NS140303BS30	0.014 inch	185cm	3.0gf	3.0cm	3cm	Straight
NS140303BJ30	0.014 inch	185cm	3.0gf	3.0cm	3cm	J
NS140303FS30	0.014 inch	300cm	3.0gf	3.0cm	3cm	Straight
NS140303FJ30	0.014 inch	300cm	3.0gf	3.0cm	3cm	J
NS180303BS25	0.018 inch	185cm	2.5gf	3.0cm	3cm	Straight
NS180303BJ25	0.018 inch	185cm	2.5gf	3.0cm	3cm	J
NS180303FS25	0.018 inch	300cm	2.5gf	3.0cm	3cm	Straight
NS180303FJ25	0.018 inch	300cm	2.5gf	3.0cm	3cm	J
NS180303BS30	0.018 inch	185cm	3.0gf	3.0cm	3cm	Straight
NS180303BJ30	0.018 inch	185cm	3.0gf	3.0cm	3cm	J
NS180303FS30	0.018 inch	300cm	3.0gf	3.0cm	3cm	Straight
NS180303FJ30	0.018 inch	300cm	3.0gf	3.0cm	3cm	J



**Add**

No.925 Jinyuan yi Road,  
Shanghai, P.R.C.



**TEL**

021-59140046



**FAX**

021-69000013



**P.C.**

201803



**WEB**

www.kdl-int.com



**Shanghai Kindly Medical Instruments Co., Ltd**