


# Pressure Bandage



- Unique design prevent unintended movement of the bandage
- Smooth transition between sheath radiopacity

Ref.NO.	Product Description	Model
FPB01	Inflatable femoral bandage, 29cm*6.5cm with 60ml syringe.	PB06
FPB01B	Inflatable femoral bandage, 31cm*6.5cm with 60ml syringe.	PB06A

# Pressure Bandage



- The transparency ensures complete visualization of the puncture site
- Dual-balloons enhance precise compression over the radial artery
- Standard luer lock connector

Ref.NO.	Product Description	Model
PB04A	Bandage:24cm*4cm Ballon:3.5cm*3.5cm	ZDX-D-18
PB04B	Bandage:27cm*4cm Ballon:3.5cm*3.5cm	ZDX-S-18

# Pressure Bandage



- Good flexibility, favorable contact
- No effects on venous blood circulation
- Pressure indication, convenient to adjust compression pressure
- Curved surface silicone available, much more comfortable for patient

Ref.NO.	Product Description	Model
PB01	Flat surface, Bandage 18cm*3cm(±0.5cm), Pressure indication 0-20N	ZXDI-18
PB02	Flat surface, Bandage 22cm*3cm(±0.5cm), Pressure indication 0-20N	ZXDI-22
PB03	Flat surface, Bandage 26cm*3cm(±0.5cm), Pressure indication 0-20N	ZXDI-26
PB01A	Curved surface, Bandage 18cm*3cm(±0.5cm), Pressure indication 0-20N	ZXDII-18
PB02A	Curved surface, Bandage 22cm*3cm(±0.5cm), Pressure indication 0-20N	ZXDII-22
PB03A	Curved surface, Bandage 26cm*3cm(±0.5cm), Pressure indication 0-20N	ZXDII-26