


Hydrophilic Guide Wire



- High quality Niti core wire enables 1:1 torque control for better navigation and faster vessel selection
- TPU jacket with tungsten providing excellent visibility under X-ray.
- Unique coating technology ensure a lubricant surface durable

Ref.NO.	Product Description		
	O.D. (inch/mm)	Length (cm)	Tip Form
HG353AF	0.035" (0.89mm)	150	45° Angled
HG354AF		180	
HG355AF		260	
HG353JF		150	J tip
HG354JF		180	
HG355JF		260	
HG383AF	0.038" (0.97mm)	150	45° Angled
HG384AF		180	
HG385AF		260	
HG383JF		150	J tip
HG384JF		180	
HG385JF		260	