

AlCath Flux G eXtra

3.5 mm Irrigated GoldTip Ablation Catheter



Ordering Information

Model	Tip	Curve	Spacing	Shaft Diameter	Length	Order Number
AlCath Flux Red G eXtra	3.5 mm irrigated	Red	2-5-2 mm	7 F	110 cm	376379
AlCath Flux Blue G eXtra	3.5 mm irrigated	Blue	2-5-2 mm	7 F	110 cm	376388
AlCath Flux Green G eXtra	3.5 mm irrigated	Green	2-5-2 mm	7 F	110 cm	376384
AlCath Flux Black G eXtra	3.5 mm irrigated	Black	2-5-2 mm	7 F	110 cm	376390
AlCath Flux Cyan G eXtra	3.5 mm irrigated	Cyan	2-5-2 mm	7 F	110 cm	376386

Technical Data

Catheter Specification	
Tip electrode length	3.5 mm
Tip electrode type	Irrigated
Tip temperature sensor	Thermocouple
Ring electrode length	1.55 mm
Number of electrodes	4
Electrode spacing	2-5-2 mm
Electrode material	Gold (Au)
Shaft diameter	7 F (2.3 mm)
Curve shape	Steerable – FullCircle
Effective catheter length	110 cm
Connection plug	9-pole connector (REDEL), female, 60° coded Luer-lock (female)
Patient cable	PK-111, 112, 142, 147, 150

Patient Cable	Ablation System	Order Number
PK-111	Osypka HAT 300 Smart	330080
PK-112	Medtronic Atakr, Medtronic Atakr II	330081
PK-142	Abbott (SJM) IBI-1500 T7 & T11 Abbott Ampere	362442
PK-147	BIOTRONIK Qubic RF, Stockert EP-Shuttle	398853
PK-150	Stockert SmartAblate	402668

Curve Type	Reach R
Red	48 mm
Blue	57 mm
Green	65 mm
Black	73 mm
Cyan	80 mm


Transport and Storage Conditions	
Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight
Shelf life	24 month



AlCath G FullCircle

4 mm GoldTip

Ablation Catheter



Ordering Information

Model	Tip	Curve	Spacing	Shaft Diameter	Length	Order Number
AlCath Red TC G FullCircle	4 mm	Red	2-5-2 mm	7 F	110 cm	364537
AlCath Blue TC G FullCircle	4 mm	Blue	2-5-2 mm	7 F	110 cm	364525
AlCath Green TC G FullCircle	4 mm	Green	2-5-2 mm	7 F	110 cm	364533
AlCath Black TC G FullCircle	4 mm	Black	2-5-2 mm	7 F	110 cm	364522
AlCath Cyan TC G FullCircle	4 mm	Cyan	2-5-2 mm	7 F	110 cm	364529

Technical Data

Catheter Specification	
Tip electrode length	4 mm
Tip electrode type	Non-irrigated
Tip temperature sensor	Thermocouple
Ring electrode length	1.55 mm
Number of electrodes	4
Electrode spacing	2-5-2 mm
Electrode material	Gold (Au)
Shaft diameter	7 F (2.3 mm)
Curve shape	Steerable - FullCircle
Effective catheter length	110 cm
Connection plug	9-pole connector (REDEL), female, 60° coded
Patient cable	PK-111, 112, 142, 147, 150

Patient Cable	Ablation System	Order Number
PK-111	Osypka HAT 300 Smart	330080
PK-112	Medtronic Atakr Medtronic Atakr II	330081
PK-142	Abbott (SJM) IBI-1500 T7 & T11 Abbott Ampere	362442
PK-147	BIOTRONIK Qubic RF, Stockert EP-Shuttle	398853
PK-150	Stockert SmartAblate	402668

Curve Type	Reach R
Red	48 mm
Blue	57 mm
Green	65 mm
Black	73 mm
Cyan	80 mm

Transport and Storage Conditions	
Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight
Shelf life	24 month

