


Lokum Amplatz Type Guide Wire

combining positioning and extra support in one design



Lokum Amplatz Type Guide Wire

- Especially designed for strong positioning and support in interventional procedures.
- High performance PTFE coating providing low surface friction and enhanced trackability.
- Intended to use for interventional peripheral and interventional biliary drainage procedures and catheter exchanges.



Technical Information

Guide Wire Material	Tempered Stainless Steel with round inner core
Coating	High performance PTFE
Diameter	0.035" (0.89 mm)
Distal Tip	Straight / 3 mm J-Tip / 6 mm J-Tip
Tip Configuration	5 cm Flexible Tip
Shaft Length	180 cm / 260 cm / 320 cm / 360 cm

* a choice of various sizes will assist your patients needs

Order Codes

A choice of various diameters and lengths.

Ø	↔	⊖	⊖	3 mm ⊖	6 mm ⊖
0.035" (0.89 mm)	180 cm	5 cm	LOAM-180-ST	LOAM-180-3T	LOAM-180-6T
0.035" (0.89 mm)	260 cm	5 cm	LOAM-260-ST	LOAM-260-3T	LOAM-260-6T
0.035" (0.89 mm)	320 cm	5 cm	LOAM-320-ST	LOAM-320-3T	LOAM-320-6T
0.035" (0.89 mm)	360 cm	5 cm	LOAM-360-ST	LOAM-360-3T	LOAM-360-6T



Andratec GmbH
 Simmener Straße. 70, D-56075 Koblenz, Germany
Office & Warehouse
 Im Deutschen Kaiser,
 Kastorstraße 3, 56068 Koblenz, Germany
 ☎ +49(0) 261 100 49 653
 📠 +49(0) 261 100 49 654
 📧 info@andratec.com
 🌐 www.andratec.com




Lokum L-Quest Type Guide Wire

balanced combination of extra strong support, high pushability and variable tip lengths



Lokum L-Quest Type Guide Wire

- Especially designed for extra strong support high pushability and catheter / delivery system exchange.
- High performance PTFE coating providing low surface friction and enhanced trackability.
- Intended to use in complex diagnostic and interventional procedures.
- 3 tips flexibilities to ensure secure and safe positioning: 5 cm, 10 cm and 15 cm.



Technical Information

Lokum L-Quest Type Guide Wire

Guide Wire Material	Tempered Tempered Stainless Steel
Coating	High performance PTFE
Diameter	0.035" (0.89 mm)
Distal Tip	Straight / 3 mm J-Tip / 6 mm J-Tip
Tip Configuration	5 cm / 10 cm / 15 cm Flexible Tip
Shaft Length	180 cm / 260 cm

Lokum L-Quest Type Mini Guide Wire

Guide Wire Material	Extra Tempered Tempered Stainless Steel
Coating	High Performance PTFE
Diameter	0.025" (0.635 mm)
Distal Tip	Straight / 3 mm J-Tip / 6 mm J-Tip
Tip Configuration	5 cm / 10 cm / 15 cm Flexible Tip
Shaft Length	150 cm / 260 cm

*a choice of various sizes will assist your patients needs

Order Codes

Lokum L-Quest Type Guide 0.035"

Ø	↔	⊖	⊕	3 mm ⊕	6 mm ⊕
0.035" (0.89 mm)	180 cm	5 cm	LOLU-05-180-ST	LOLU-05-180-3T	LOLU-05-180-6T
0.035" (0.89 mm)	180 cm	10 cm	LOLU-10-180-ST	LOLU-10-180-3T	LOLU-10-180-6T
0.035" (0.89 mm)	180 cm	15 cm	LOLU-15-180-ST	LOLU-15-180-3T	LOLU-15-180-6T
0.035" (0.89 mm)	260 cm	5 cm	LOLU-05-260-ST	LOLU-05-260-3T	LOLU-05-260-6T
0.035" (0.89 mm)	260 cm	10 cm	LOLU-10-260-ST	LOLU-10-260-3T	LOLU-10-260-6T
0.035" (0.89 mm)	260 cm	15 cm	LOLU-15-260-ST	LOLU-15-260-3T	LOLU-15-260-6T

Lokum L-Quest Mini Type Guide 0.025"

Ø	↔	⊖	⊕	3 mm ⊕	6 mm ⊕
0.025" (0.635 mm)	150 cm	5 cm	LUMI-05-150-ST	LUMI-05-150-3T	LUMI-05-150-6T
0.025" (0.635 mm)	150 cm	10 cm	LUMI-10-150-ST	LUMI-10-150-3T	LUMI-10-150-6T
0.025" (0.635 mm)	150 cm	15 cm	LUMI-15-150-ST	LUMI-15-150-3T	LUMI-15-150-6T
0.025" (0.635 mm)	260 cm	5 cm	LUMI-05-260-ST	LUMI-05-260-3T	LUMI-05-260-6T
0.025" (0.635 mm)	260 cm	10 cm	LUMI-10-260-ST	LUMI-10-260-3T	LUMI-10-260-6T
0.025" (0.635 mm)	260 cm	15 cm	LUMI-15-260-ST	LUMI-15-260-3T	LUMI-15-260-6T



Andratec GmbH
 Simmener Straße. 70, D-56075 Koblenz, Germany
Office & Warehouse
 Im Deutschen Kaiser,
 Kastorstraße 3, 56068 Koblenz, Germany
 ☎ +49(0) 261 100 49 653
 📠 +49(0) 261 100 49 654
 📧 info@andratec.com
 🌐 www.andratec.com

