

MultiCath 4-Pole

Multipolar Fixed Curve Catheter for Electrophysiological Diagnostics



Ordering Information


Model	Curve	Spacing	Shaft Diameter	Length	Order Number
MultiCath 4J Sp ø4F 5 mm	Josephson special	5 mm	4 F	100 cm	356728
MultiCath 4J Sp ø4F 2-5-2 mm	Josephson special	2-5-2 mm	4 F	100 cm	356727
MultiCath 4J Sp ø5F 2 mm	Josephson special	2 mm	5 F	80 cm	351170
MultiCath 4J Sp 100 cm ø5F 5 mm	Josephson special	5 mm	5 F	100 cm	356434
MultiCath 4J Sp ø5F 2-5-2 mm	Josephson special	2-5-2 mm	5 F	80 cm	351171
MultiCath 4J Sp 100 cm ø5F 2-5-2 mm	Josephson special	2-5-2 mm	5 F	100 cm	356433
MultiCath 4C ø5F 2-5-2 mm	Cournand	2-5-2 mm	5 F	80 cm	351177
MultiCath 4HIS ø5F 5 mm	HIS	5 mm	5 F	100 cm	359385
MultiCath 4J ø6F 5 mm	Josephson	5 mm	6 F	110 cm	351185
MultiCath 4J ø6F 5 mm SoftTip	Josephson	5 mm	6 F	110 cm	351182
MultiCath 4D ø6F 5 mm	Damato	5 mm	6 F	110 cm	351189
MultiCath 4C ø6F 5 mm	Cournand	5 mm	6 F	110 cm	351188
MultiCath 4C ø6F 5 mm SoftTip	Cournand	5 mm	6 F	110 cm	351183
MultiCath 4M ø6F 5 mm	Multipurpose	5 mm	6 F	110 cm	351194

MultiCath 4-Pole

Technical Data

Catheter Specification		
Tip electrode length	2.5 mm (4 F & 5 F)	3.2 mm (6 F & 6 F Soft Tip)
Tip electrode type	Standard	
Ring electrode length	1.5 mm (4 F & 5 F)	2.2 mm (6 F) 1.3 mm (6 F Soft Tip)
Number of electrodes	4	
Electrode spacing	2 mm, 5 mm, 2-5-2 mm	
Electrode material	Platinum Iridium (Pt/Ir)	
Shaft diameter	4F (1.3 mm), 5F (1.7 mm), 6F (2.0 mm)	
Curve type	Josephson, Josephson special, Cournand, Damato, HIS, Multipurpose	
Curve shape	Fixed (Standard or Soft Tip)	
Effective catheter length	80 cm, 100 cm, 110 cm	
Connection plug	4-pole connector (REDEL), female, Encoding 0°	
Patient cable	MPK-4-R	

Transport and Storage Conditions	
Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight
Shelf life	24 months



MultiCath 10-Pole

Multipolar Fixed Curve Catheter for Electrophysiological Diagnostics



Ordering Information

Model	Curve	Spacing	Shaft Diameter	Length	Order Number
MultiCath 10J Sp ø5F 2-5-2 mm	Josephson special	2-5-2 mm	5 F	100 cm	359388
MultiCath 10J Sp ø5F 2-8-2 mm	Josephson special	2-8-2 mm	5 F	100 cm	359389
MultiCath 10J ø6F 5 mm	Josephson	5 mm	6 F	110 cm	351193
MultiCath 10J CS ø6F 2-5-2 mm	Josephson	2-5-2 mm	6 F	110 cm	351192
MultiCath 10J CS ø6F 2-8-2 mm	Josephson	2-8-2 mm	6 F	110 cm	351191

Technical Data

Catheter Specification

Tip electrode length	2.5 mm (5 F)	3.2 mm (6 F Soft Tip)
Tip electrode type	Standard	
Ring electrode length	1.5 mm (5 F)	1.3 mm (6 F Soft Tip)
Number of electrodes	10	
Electrode spacing	2-5-2 mm, 2-8-2 mm, 5 mm	
Electrode material	Platinum Iridium (Pt/Ir)	
Shaft diameter	5 F (1.7 mm), 6 F (2.0 mm)	
Curve type	Josephson, Josephson special	
Curve shape	Fixed (Standard or Soft Tip)	
Effective catheter length	100 cm, 110 cm	
Connection plug	10-pole connector (REDEL), female, Encoding 0°	
Patient cable	MPK-10-R	

Transport and Storage Conditions

Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight
Shelf life	24 months

MultiCath 10
Josephson

MultiCath 10
Josephson Special

