




**PROLONG™**


**Guide Extension Catheter**



Where there is a **Core**  
There is a way to  
**Guide**

## Two Radiopaque Marks

Reinforce visibility and locating



### Soft Tip


Minimizes damage to vessels

### Stainless Hypo Tube

Increases pushability and folding resistance, and reduced device interaction

### Hydrophilic Coating

Improves conveyability in complex lesions




### Helical Collar

Maintains overall solidity and strength



### Unique Lining Core Design

Enhances supporting performance



#### Ordering Information

Ref. No.	Size	Guide Extension Catheter				Lining Core		
		Working Length	Tubing Length	I.D.	O.D.	Lined Tube Length	I.D.	O.D.
GEC015F25	5.0Fr	150cm	25cm	1.17mm	1.35mm	28mm	0.40mm	1.12mm
GEC015F40	5.0Fr	150cm	40cm	1.17mm	1.35mm	43mm	0.40mm	1.12mm
GEC0155F25	5.5Fr	150cm	25cm	1.30mm	1.60mm	28mm	0.40mm	1.25mm
GEC0155F40	5.5Fr	150cm	40cm	1.30mm	1.60mm	43mm	0.40mm	1.25mm
GEC016F25	6.0Fr	150cm	25cm	1.45mm	1.70mm	28mm	0.40mm	1.40mm
GEC016F40	6.0Fr	150cm	40cm	1.45mm	1.70mm	43mm	0.40mm	1.40mm
GEC017F25	7.0Fr	150cm	25cm	1.57mm	1.90mm	28mm	0.40mm	1.52mm
GEC017F40	7.0Fr	150cm	40cm	1.57mm	1.90mm	43mm	0.40mm	1.52mm
GEC018F25	8.0Fr	150cm	25cm	1.80mm	2.16mm	28mm	0.40mm	1.75mm
GEC018F40	8.0Fr	150cm	40cm	1.80mm	2.16mm	43mm	0.40mm	1.75mm
GEC025F25	5.0Fr	150cm	25cm	1.17mm	1.35mm	28mm	0.40mm	1.12mm
GEC025F40	5.0Fr	150cm	40cm	1.17mm	1.35mm	43mm	0.40mm	1.12mm
GEC0255F25	5.5Fr	150cm	25cm	1.30mm	1.60mm	28mm	0.40mm	1.25mm
GEC0255F40	5.5Fr	150cm	40cm	1.30mm	1.60mm	43mm	0.40mm	1.25mm
GEC026F25	6.0Fr	150cm	25cm	1.45mm	1.70mm	28mm	0.40mm	1.40mm
GEC026F40	6.0Fr	150cm	40cm	1.45mm	1.70mm	43mm	0.40mm	1.40mm
GEC027F25	7.0Fr	150cm	25cm	1.57mm	1.90mm	28mm	0.40mm	1.52mm
GEC027F40	7.0Fr	150cm	40cm	1.57mm	1.90mm	43mm	0.40mm	1.52mm
GEC028F25	8.0Fr	150cm	25cm	1.80mm	2.16mm	28mm	0.40mm	1.75mm
GEC028F40	8.0Fr	150cm	40cm	1.80mm	2.16mm	43mm	0.40mm	1.75mm

#### Shanghai Kindly Medical Instruments Co.,Ltd.

ADD: No.925, Jinyuan Yi Road, Shanghai, China

TEL: +86-21-59140056

P.C.: 201803

FAX: +86-21-69000013

https://sh-intmedical.com/