

ViaCath 10-Pole

Steerable Diagnostic Catheter



Ordering Information

Model	Curve	Spacing	Shaft Diameter	Length	Order Number
ViaCath NG 10/S/2-6-2 mm	Standard	2-6-2 mm	6 F	110 cm	351200
ViaCath NG 10/L/2-6-2 mm	Large	2-6-2 mm	6 F	110 cm	358797
ViaCath NG 10/L/2-8-2 mm	Large	2-8-2 mm	6 F	110 cm	370309
ViaCath NG 10/XL/2-10-2 mm	Extra Large	2-10-2 mm	6 F	110 cm	370144


Technical Data

Catheter Specification	
Tip electrode length	2 mm
Tip electrode type	Standard
Ring electrode length	1 mm
Number of electrodes	10
Electrode spacing	2-6-2, 2-8-2, 2-10-2 mm
Electrode material	Platinum Iridium (Pt/Ir)
Shaft diameter	6 F (2.0 mm)
Curve type	Standard, Large, Extra Large
Curve shape	Steerable
Effective catheter length	110 cm
Connection plug	10-pole connector (REDEL), female, 0° coded
Patient cable	MPK-10-R

Transport and Storage Conditions	
Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight
Shelf life	24 months

ViaCath 20-Pole

Steerable Diagnostic Catheter



Ordering Information

Model	Curve	Spacing	Shaft Diameter	Length	Order Number
ViaCath 20/XL/2-8-2 mm	Extra Large	2-8-2 mm	6 F	90 cm	351201

Technical Data

Catheter Specification	
Tip electrode length	2 mm
Tip electrode type	Standard
Ring electrode length	1 mm
Number of electrodes	20
Electrode spacing	2-8-2 mm
Electrode material	Platinum Iridium (Pt/Ir)
Shaft diameter	6 F (2.0 mm)
Curve type	Extra Large
Curve shape	Steerable
Effective catheter length	90 cm
Connection plug	2 x 10-pole connector (REDEL), female, 0° coded
Patient cable	2 x MPK-10-R

Transport and Storage Conditions	
Storage	Please store this sterile product in the original sterile packaging
Light	Keep away from sunlight
Shelf life	24 months